

HIOSSSEN IS3™ ISQ

HIOSSSEN



Rev. 2018/07/09

MARKETING DEPT.

PM

Contents

- I. Introduction
- II. Special Feature
- III. Practical Use
- IV. COMPETITOR COMPARISON

I. Introduction



Diagnostic & Examination Tools



Visual,
Percussion



Imaging




Magnetic Resonance Measuring
Device - ISQ

Interproximal dental caries



Percussion

“ The method of tapping the patient’s body to check for abnormalities ”



Percussion

“Percussion” is used as a method to diagnosis the abnormality of a knee, internal organs or teeth, but a more specialized method should be used for accurate diagnosis.




ISQ ; Impant Stability Quotient

HIOSSEN IS3




ISQ Measurement methodology



The technique is Resonance Frequency Analysis (RFA) and the measurement unit is Implant Stability Quotient (ISQ).

The Multipeg™ is excited by magnetic pulses and vibrates through the stiffness in the contact area between the bone and the implant surface.

High ISQ Value = Reduce Risk / Reduce Healing time




Values in the range of 75 ISQ means the implant is already so stable that osseointegration cannot add stability. The proof of osseointegration is the lack of a drop in ISQ.

If the implant has a low initial ISQ, let's say 55, osseointegration will add stability over time.

ISQ correlates to micro mobility.

The ISQ-scale 1-99 ISQ

It seems that ISQ-values above 70 can withstand normal forces in the mouth.



Values above 70 ISQ indicates a very stable implant with low micro mobility . Typically used for one-stage and immediate loading.

A second measurement, before the final restoration is recommended to verify osseointegration.

How the stability is varied

Each implant is moulded into plaster and then placed into plastic cylinders with only the top of the implant showing.


The cylinders are mounted into a device that can vary the stability through the whole ISQ-scale in a consistent, predictable, repeatable way.



II. IS3 Special Feature

1. Treatment Guidance & Predictability

Assists clinicians in the treatment of patients with risk factors



*Ref) Loma linda. Immediate loading, JOI, Vol. XXXVIII, Special Issue no. One/2012

2. Ease of use & Objective, non-invasive operation

By mounting a Multipeg™ the measurement is made in a second.

Just aim for the magnet on top of the Multipeg™ (as shown below)

Non-invasive, objective, accurate and repeatable.



3. Reusable MultiPeg™ with various implant brand platform


Determines optimal loading protocol

- Tissue friendly, durable titanium
- Sealed magnet
- Autoclavable app. 20 times
- Laser marked with type number
- Optimal platform fit
- ISQ Standard Calibrated™



4. Optimized Patient Care

Determines optimal loading protocol



5. High Accuracy measurement results

ISQ Standard Calibrated™ (ISC)




-Optimal physical fit

Any misfit in the platform connection will be exposed by comparing against ISC.

-Minimum variance

By elaborating the physical dimensions of the Multipeg™ it is possible to minimize the variance against the ISC.




Empirical RFA data, mainly from Smartpegs type 1 and 4, was used to create an ISQ Standard Calibration graph as a test reference for Multipeg™.

5. High Accuracy measurement results




MulTipeg™ (Type #74)

VS

Smartpeg (Type #50)



This "border" is supposed to fit on the top of the implant platform



Prolonged Shank + gap

6. Scientifically verified



More than 700 articles on the subject have been published in peer-reviewed papers since 1996.

**Resonance Frequency Analysis for Implant Stability Measurements,
INTEGRATION DIAGNOSTICS UPDATE 2015;1:1-11*

III. IS3 Practical Use

1) Package Components (Purchase Includes)



IS3™ Unit

- *Recommended initial battery charging time : 3 hours
- *2Years manufacturer warranty



MulTipeg™ Driver

- *Guaranteed for at least 100 autoclave cycle



MulTipeg™

- *ET Standard connection included
- *More option
Refer to the assortment to find the compatible implant platform



Main adapter & insert

- *Compatibility with worldwide power sockets

2) Measurement



STEP 1

Select the correct **Multipeg™** for the implant being measured.

STEP 2

Attach the Multipeg™ using the driver included. **Hand tighten (6~8Ncm).**



STEP 3

Turn on the instrument and hold the tip of the instrument close to the top of the Multipeg™.

STEP 4

The instrument emits a beeping sound when it starts measuring and another sound when the ISQ value is presented on the display a second later.



3) Cleaning

- **IS3™ Unit**

The unit is cleaned with a wet cloth. Any disinfectants can be used.

- **MulTipeg™ and MulTipeg™ Driver**

It should be cleaned with water or detergent using a light brush. For use in environments requiring sterility, the MulTipeg™ and driver should be autoclaved before use according to below instructions.

At least 3 minutes at 134°C (273°F)



4) Troubleshooting

- **Difficult to achieve a measurement**



Monitor closely – should not touch soft tissue

- **Noise warning (audible & visible on the display)**



An electric device close to the unit is causing the warning. Try to isolate and remove the source.

- **Unit turns off suddenly**

- Automatically turn off after 1min of inactivity
- Battery level is too low.

- **Error code (if malfunctioning, these codes are shown on the display)**

E1 Hardware error. Malfunctioning electronics

E2 Noise error. Shown if constant electromagnetic noise is present

E3 Pulse power error. Malfunctioning magnetic pulse generation

5) Indications

- 1) Right after implantation
- 2) 2nd surgery
- 3) Every 3weeks to 1month during osseointegration
- 4) When prosthesis prepared
- 5) Final prosthesis applied

* By monitoring repeatedly of the patient's ISQ value, optimal loading timing can be determined. If the ISQ value is 70 or higher, it can be considered as prosthetic loading time.

IV. COMPETITOR COMPARISON

1. Competitor comparison

	HIOSSSEN IS3	OSSTELL ISQ	
Image			
Design	 Portable pen type	Probe with Docking station	Light, Mobility
Probe Material	 PEEK	Stainless steel	Biocompatible material
Measurement Technique	Resonance Frequency Analysis	Resonance Frequency Analysis	-
Peg part	 Reusable(up to 20times)	Single use	More saving, Durability

The End.