OSSTEM[®]
IMPLANT

123 KIT

User Manual

OSSTEM[®]
IMPLANT

CE
2460

Storage Condition :
Dry place at room temperature(1°C~30°C)

Caution, Consult
accompanying
Documents

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English

■ Features of the 123 KIT

1. The 123 Kit has the drilling sequences dramatically reduced for use in conventional fixture installation.
2. The 123 Kit has drills with excellent cutting capabilities to minimize the heat generated from drilling.
3. This Kit also has a drill stop which will give excellent control of the drill depth.

■ About the Color Coding

1. The color coding of the drills are determined by the diameter of the fixture, as follows: 3.5 Fixtures (Yellow), 4.0 Fixtures (Green), 4.5 Fixtures (Blue) and 5.0 Fixtures (Red).
2. The Color coding for the drill stops are classified by the length: 6mm (Purple), 7mm (Gray), 8mm (Yellow), 9.5mm (Red), 11mm (Blue), 12.5mm (Green), 14mm (Black) and 16mm (Purple).


■ Cleaning the KIT

Clean the Surgical tool after each use as follows:

1. After surgery, clean the tools using the conventional method, i. e. by soaking them in alcohol.
2. After cleaning the tools, remove any remaining blood stains or foreign matter completely with distilled water or under running water.
3. Remove moisture completely using a dry cloth or a heat fan.
4. Insert and secure the moisture free surgical tools into the correct slots of the kit.
5. After placing the tools in the kit, sterilize in the autoclave (132°C for 15 minutes) and store at room temperature.

■ Cautions

1. When drilling, move the hand piece perpendicularly up and down in a pumping motion.
2. To reduce the friction during drilling, use a constant spray of saline solution.

[Figure. Baseline of Laser Marking]

3. The laser marking indicates the drilling depth levels. The bottomline of the markings serves as the baseline for depth. On all tools the laser marking for 10mm and 11.5mm are marked with a thick line.
4. The drill stop should be used because the cutting force of the 123 drill is so high it is hard to control the depth.
5. The actual length of the drill stop is defined after connecting to the drill.
6. The drill can be used for up to 50 times.
7. The maximum allowable torque of the 0.9 hex driver is 20 Ncm. And that of the 1.2 hex driver is 35 Ncm. The driver can also be used up to 50 times.
8. To connect the drill to the drill extension, let the flat handle tip face the opposite side of the slot. Avoid excessive force during the connection.
9. Separate all used tools immediately after surgery and after cleaning and drying, store at room temperature.
10. Do not leave the instruments in a place vulnerable to contamination.
11. Be sure to sterilize surgical tools in the autoclave (132°C for 15 minutes) before each new surgical procedure.
12. Do not use hydrogen peroxide as a disinfectant or cleaner (it may cause damage or discoloration of the Tin coating, laser marking and color coding).
13. This product comes with a 1 year warranty from the date of purchase.

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OBELIS S.A.