


Drilling Sequence **Straight Drill**

TSIII Ultrawide CA/SA | **SSIII Ultrawide CA/SA** | **USIII Ultrawide CA/SA** (Length : 10mm)

Ø6.0mm



Bone Quality	Ø 2.2 Drill	Ø 2.0/3.0 Pilot Drill	Ø 3.0 Drill	Ø 3.0/3.8 Pilot Drill	Ø 3.8 Drill	Ø 4.6 Drill	Ø 5.2 Direct drill	Ø 5.5 Direct drill	F6.0 Cortical Drill	Ø 6.0 Fixture
Soft	▶	▶	▶	▶	▶	▶	▶			
Normal	▶	▶	▶	▶	▶	▶	▶	▶		Implant placement
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	

Ø7.0mm




Bone Quality	Ø 2.2 Drill	Ø 2.0/3.0 Pilot Drill	Ø 3.0 Drill	Ø 3.0/3.8 Pilot Drill	Ø 3.8 Drill	Ø 4.6 Drill	Ø 5.5 Direct drill	Ø 6.2 Direct drill	Ø 6.5 Direct drill	F7.0 Cortical Drill	Ø 7.0 Fixture
Soft	▶	▶	▶	▶	▶	▶	▶	▶			
Normal	▶	▶	▶	▶	▶	▶	▶	▶	▶		Implant placement
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	▶	

Drilling Sequence **Taper Drill**


TSIII Ultrawide CA/SA | **SSIII Ultrawide CA/SA** | **USIII Ultrawide CA/SA** (Length : 10mm)

Ø6.0mm



Bone Quality	Ø 2.2 Drill	Ø 2.0/3.0 Pilot Drill	Ø 3.0 Drill	Ø 3.0/3.8 Pilot Drill	Ø 3.8 Drill	F6.0 Taper Drill	F6.0 Taper Drill	F6.0 Cortical Drill	Ø 6.0 Fixture
Soft	▶	▶	▶	▶		▶			
Normal	▶	▶	▶	▶	▶		▶		Implant placement
Hard	▶	▶	▶	▶	▶	▶	▶	▶	

Ø7.0mm



Bone Quality	Ø 2.2 Drill	Ø 2.0/3.0 Pilot Drill	Ø 3.0 Drill	Ø 3.0/3.8 Pilot Drill	Ø 3.8 Drill	F6.0 Taper Drill	F7.0 Taper Drill	F7.0 Taper Drill	F7.0 Cortical Drill	Ø 7.0 Fixture
Soft	▶	▶	▶	▶		▶	▶			
Normal	▶	▶	▶	▶	▶	▶		▶		Implant placement
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	

Recommended insertion torque : ≤40Ncm

※Recommend reducing the insertion speed to 15rpm or lower, due to the TSIV/USIV aggressive threads.

Recommended insertion torque: ≤40Ncm

TS fixture insertion depth: In normal/hard bone, the implant should be placed 1mm subcrestal. In soft bone, it should be inserted at bone level for higher stability.