


ZioCera

Abutment



Feature

- Cement / screw / combination retained prosthesis
- Single / bridge restoration / anterior area
(Not recommended : posterior area case)
- Fixture level impression
- Zirconia material appropriate for all ceramic prosthesis fabrication in anterior area
- Natural dentin color abutment shade
- Bio friendly and excellent strength
- 2 types : better surgery convenience (straight / 17° angled)
- Use exclusive abutment screw
- Material : zirconia (non coating) / Ti-6Al-4V (WCC coating)
- Connect using 1.2 hex torque driver
- Recommended tightening torque : mini 20Ncm / regular 30Ncm

Straight


D	Ø 4.5 / 5.5 / 6.5 mm
G/H	3.5 / 5.0 mm
Type	Hex / Non-Hex

Angled

D	Ø 5.5 / 6.5 mm
G/H	3.0 / 4.0 mm
Type	Hex / Non-Hex


Cement Retained Restoration

- Cement retained type of all ceramic prosthesis fabrication is recommended for ZioCera abutment
- Use zirconium exclusive bur for modifying abutment, and use irrigation




Screw Retained Restoration

- Screw retained type prosthesis fabrication is possible with direct build up
- Use zirconium exclusive bur for modifying abutment, and use irrigation
- Fabrication of esthetic implant prosthesis is possible with exclusive porcelain powder build up



Prosthetic Flow Diagram



Prosthetic Process

Fixture Level Impression Cement Type Prosthesis

01

Remove healing abutment

- Remove healing abutment with 1.2 hex hand driver by hand



1.2 Hex Hand Driver



02

Impression

- Consider abutment diameter and type (hex/non-hex)
- Select impression coping specification (pick-up / transfer type)
- Connect using 1.2 hex hand driver by hand
- Block out driver hole of transfer impression coping
- Perform peri apical x-ray to check correct connection
- Take impression by applying impression material around impression coping first




Fixture Transfer Impression Coping



03 Lab Side

Fabricate working model

- Check impression body
- Fabricate working model in normal way and pour stone
- Connect Ziocera abutment and modify path and customize shape
- Use Zirconia exclusive bur
- Must use irrigation while modifying (High heat generated while modifying can break abutment)



04 Lab Side

Fabricate ceramic coping

- Fabricate ceramic coping in normal way



05 Lab Side

Porcelain build up

- Porcelain build up and firing on ceramic coping
- Polishing procedure in normal way
- Check prosthesis in the working model



06 Lab Side

Abutment connection

- Check delivered prosthesis from the lab
- Remove healing abutment or temporary prosthesis from mouth
- Connect using 1.2 hex torque driver (mini 20Ncm / regular 30Ncm)
- Check right connection with x-ray



1.2 Hex Torque Driver



Torque Wrench



07 Lab Side

Connect final prosthesis

- Abutment screw hole block out
- Connect prosthesis by cementation and remove cement completely



Fixture Level Impression Screw Type Prosthesis

01

Remove healing abutment

- Remove healing abutment with 1.2 hex hand driver by hand



1.2 Hex Hand Driver



02

Impression

- Consider abutment diameter and type (hex/non-hex)
- Select impression coping specification (pick-up / transfer type)
- Connect using 1.2 hex hand driver by hand
- Block out driver hole of transfer impression coping
- Perform peri apical x-ray to check correct connection
- Take impression by applying impression material around impression coping first




1.2 Hex Hand Driver



03 Lab Side

Fabricate working model


- Check impression body
- Fabricate working model in normal way and pour stone
- Connect Ziocera abutment and modify path and customize shape
- Use Zirconia exclusive bur
- Must use irrigation while modifying (High heat generated while modifying can break abutment)



04 Lab Side

Porcelain build up and firing

- Porcelain build up with Zirconia exclusive powder
- Easy to form screw hole using waxing screw for lab
- To prevent change of mechanical property, limit firing to 5 times
- Polishing procedure in normal way
- Check prosthesis in the working model



05

Connect final prosthesis

- Check delivered prosthesis from the lab
- Remove healing abutment or temporary prosthesis from mouth
- Connect using 1.2 hex torque driver (mini 20Ncm / regular 30Ncm)
- Check right connection with x-ray
- Block out Screw hole with resin



1.2 Hex Torque Driver



Torque Wrench



※ Cautions for Zirconia abutment Use

- ① Use Zirconia exclusive bur
- ② Must irrigate while milling to prevent overheating
- ③ Apply round shape to edge or corner to prevent fracture
- ④ Use zirconia exclusive power for build up