



# SURGICAL KIT



**178** 122 Taper KIT  
**179** 122 Taper Full KIT  
**186** Taper KIT  
**187** Taper Ultra KIT  
**194** 123 Straight Simple KIT  
**195** 123 Straight KIT  
**196** 123 Straight Full KIT  
**204** Ultra KIT  
**216** 485 KIT  
**220** Assist KIT  
**222** Surgical Instrument  
**244** Prosthetic Simple KIT  
**245** Prosthetic KIT  
**254** CAS KIT  
**260** LAS KIT  
**261** LAS Full KIT  
**264** Denture 4U KIT  
**270** Positioning Guide KIT  
**271** Positioning Guide Full KIT

**275** SmartGuide KIT  
**278** ESSET KIT  
**282** IM-Cure KIT  
**286** ESR KIT  
**287** ESR Full KIT  
**296** EFR KIT  
**297** EFR Full KIT  
**302** Dr.Cho's Instrument KIT  
**303** Osstem Basic Instrument KIT  
**306** Custom KIT  
**307** Healing Case  
**308** Osteo KIT  
**309** Osteotome KIT  
**310** Sinus KIT  
**311** Bone Spreader KIT  
**312** Ridge Split KIT **Straight**  
**313** Ridge Split KIT **Offset**

122 Taper KIT (O122TPK)

Applicable Products TSIII / IV SSIII



Top panel components



122 Taper Full KIT (O122TPFK)

Applicable Products TSIII / IV SSIII III / IV Ultra-wide

Top panel components



• More details on KIT components can be found in Surgical Instruments(222p~242p)


• More details on KIT components can be found in Surgical Instruments(222p~242p)

Drilling Sequence **122 Taper Drill**

**TSIII | SSIII**


(Length : 10mm)

**Ø3.0mm**




Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	122Taper Drill (F3.5)	Ø3.0 Implant
Soft	▶			Implant Placement
Normal	▶	▶		
Hard	▶		▶	

**Ø3.5mm**



Bone Quality	Sidecut Drill	122Taper Drill (F3.5)	122Taper Drill (F4.0)	Ø3.5 Implant
Soft	▶			Implant Placement
Normal	▶	▶		
Hard	▶	▶	▶	

**Ø4.0mm**




Bone Quality	Sidecut Drill	122Taper Drill (F3.5)	122Taper Drill (F4.0)	122Taper Drill (F4.5)	Ø4.0 Implant
Soft	▶	▶			Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶		▶	

F5.5 taper cortical drill marking line : Bottom line for placing 6mm or smaller implants, midline for 7mm implants, and top line for 8,5mm or greater implants

Recommended placement torque ≤ 40Ncm


TS implant placement depth : For normal bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

**Ø4.5mm**




Bone Quality	Sidecut Drill	122Taper Drill (F3.5)	122Taper Drill (F4.0)	122Taper Drill (F4.5)	122Taper Drill (F5.0)	Ø4.5 Implant
Soft	▶			▶		Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶			▶	

**Ø5.0mm**



Bone Quality	Sidecut Drill	122Taper Drill (F3.5)	122Taper Drill (F4.5)	122Taper Drill (F5.0)	122Taper Drill (F5.5)	Ø5.0 Implant
Soft	▶			▶		Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶			▶	


**Ø5.5mm**




Bone Quality	Sidecut Drill	122Taper Drill (F3.5)	122Taper Drill (F5.0)	122Taper Drill (F5.5)	Taper Cortical Drill (F5.5)	Ø5.5 Implant
Soft	▶			▶		Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶		▶	▶	

Drilling Sequence **122 Taper Drill**  
**TSIII Ultra-wide | SSIII Ultra-wide**  
(Length : 10mm)

Ø6.0mm




Ø7.0mm




F5.5 taper cortical drill marking line : Bottom line for placing 6mm or smaller implants, midline for 7mm implants, and top line for 8,5mm or greater implants  
Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

Drilling Sequence **122 Taper Drill**  
**TSIV**  
(Length : 10mm)


Ø4.0mm



Ø4.5mm




Ø5.0mm




Drilling Sequence **122 Taper Drill**  
**TSIV Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Bone Quality	Sidecut Drill	122 Taper Drill (F5.0)	122 Taper Drill (F5.5)	Ø6.0 Implant
D4	▶	▶		Implant Placement
Soft	▶	▶	▶	

Ø7.0mm



Bone Quality	Sidecut Drill	122 Taper Drill (F5.0)	122 Taper Drill (F5.5)	122 Taper Drill (F6.0)	Ø7.0 Implant
D4	▶	▶	▶		Implant Placement
Soft	▶	▶	▶	▶	

OSSTEM<sup>®</sup>  
IMPLANT

F5.5 taper cortical drill marking line : Bottom line for placing 6mm or smaller implants, midline for 7mm implants, and top line for 8,5mm or greater implants  
Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

Taper KIT (OTSK)


Applicable Products TSIII / IV SSIII

Top panel components

Torque Wrench  
TW30B



Depth Gauge  
OSDG



Taper Ultra KIT (HULTPK)


Applicable Products III Ultra-wide

Top panel components

Open Wrench  
SPOW



Ratchet Wrench  
RCWC




• More details on KIT components can be found in Surgical Instruments(222p~242p)

• More details on KIT components can be found in Surgical Instruments(222p~242p)


Drilling Sequence **Taper Drill**  
**TSIII | SSIII**  
(Length : 10mm)

Ø3.0mm




Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Cortical Drill (F3.0)	Ø3.0 Implant
Soft	▶			Implant Placement
Normal	▶	▶		
Hard	▶	▶	▶	

Ø3.5mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Cortical Drill (F3.5)	Taper Implant Tap (F3.5)	Ø3.5 Implant
Soft	▶	▶				Implant Placement
Normal	▶		▶			
Hard	▶		▶	▶		
Hard (Option)	▶		▶		▶	


Ø4.0mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.0)	Taper Cortical Drill (F4.0)	Taper Implant Tap (F4.0)	Ø4.0 Implant
Soft	▶	▶				Implant Placement
Normal	▶	▶	▶			
Hard	▶	▶	▶	▶		
Hard (Option)	▶	▶	▶		▶	


Taper cortical drill marking line : Bottom line for placing 8,5mm or greater implants, and top line for 10mm or greater implants  
Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability  
For implant tap used in hard bones, engine (25rpm recommended) is used or a torque wrench is used after assembling mount extension  
(F5.0 Implant Tap : Bottom line for placing 7.0mm or smaller implants, and top line for 8.5mm or greater implants)

Ø4.5mm




Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.0)	Taper Drill (F4.5)	Taper Cortical Drill (F4.5)	Taper Implant Tap (F4.5)	Ø4.5 Implant
Soft	▶	▶	▶				Implant Placement
Normal	▶	▶		▶			
Hard	▶	▶		▶	▶		
Hard (Option)	▶	▶		▶		▶	

Ø5.0mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.5)	Taper Drill (F5.0)	Taper Cortical Drill (F5.0)	Taper Implant Tap (F5.0)	Ø5.0 Implant
Soft	▶	▶	▶				Implant Placement
Normal	▶	▶	▶	▶			
Hard	▶	▶	▶	▶	▶		
Hard (Option)	▶	▶	▶	▶		▶	


Ø5.5mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.5)	Taper Drill (F5.5)	Taper Implant Tap (F5.5)	Ø5.5 Implant
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	


Drilling Sequence **Taper Drill**  
**TSIV**  
(Length : 10mm)

Ø4.0mm




Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Ø4.0 Implant
D4	▶		Implant Placement
Soft	▶	▶	

Ø4.5mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Drill (F4.0)	Ø4.5 Implant
D4		▶			Implant Placement
Soft	▶		▶	▶	


Ø5.0mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Drill (F4.5)	Ø5.0 Implant
D4		▶			Implant Placement
Soft	▶		▶	▶	


Drilling Sequence **Taper Drill**  
**TSIII Ultra-wide** | **SSIII Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Taper Drill (F6.0)	Taper Drill (F6.0)	Cortical Drill (F6.0)	Ø6.0 Implant
Soft	▶	▶	▶	▶		▶			Implant Placement
Normal	▶	▶	▶	▶	▶		▶		
Hard	▶	▶	▶	▶	▶		▶	▶	

Ø7.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Taper Drill (F6.0)	Taper Drill (F7.0)	Taper Drill (F7.0)	Cortical Drill (F7.0)	Ø7.0 Implant
Soft	▶	▶	▶	▶		▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		▶		
Hard	▶	▶	▶	▶	▶	▶		▶	▶	


Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal/hard bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

Drilling Sequence **Taper Drill**  
**TSIV Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Ø7.0mm





OSSTEM<sup>®</sup>  
IMPLANT

Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal/hard bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability



123 Straight Simple KIT (O123K)

Applicable Products TSII / III SSII / III



123 Straight KIT (O123FK)

Applicable Products TSII / III SSII / III



• More details on KIT components can be found in Surgical Instruments(222p~242p)

• More details on KIT components can be found in Surgical Instruments(222p~242p)

# 123 Straight Full KIT (O123STFK)


Applicable Products TSII / III SSII / III II / III Ultra-wide

Top panel components

Torque Wrench  
TW30B




Depth Gauge  
OSDG



• More details on KIT components can be found in Surgical Instruments(222p~242p)


Drilling Sequence **II Type 123 Twist Drill**  
**TSII | SSII**  
(Length : 10mm)

Ø3.5mm




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Cortical Drill	Ø3.5 Implant
Soft	▶	▶		Implant Placement
Normal	▶	▶		
Hard	▶	▶	▶	

Ø4.0mm




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Cortical Drill	Ø4.0 Implant
Soft	▶	▶			Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶	▶	▶	

Ø4.5mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Twist Drill (Ø3.0/4.1)	Cortical Drill	Ø4.5 Implant
Soft	▶		▶			Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶		▶	▶	


Ø5.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.1)	Twist Drill (Ø3.0/4.6)	Cortical Drill	Ø5.0 Implant
Soft	▶		▶			Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶		▶	▶	


Drilling Sequence **III Type 123 Twist Drill**  
**TSIII | SSIII**  
(Length : 10mm)

Ø3.5mm




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Cortical Drill (F3.5) Bottom line	Cortical Drill (F3.5) Top line	Ø3.5 Implant
Soft	▶	▶			Implant Placement
Normal	▶	▶	▶		
Hard	▶	▶		▶	

Ø4.0mm




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Cortical Drill (F4.0) Bottom line	Cortical Drill (F4.0) Top line	Ø4.0 Implant
Soft	▶	▶				Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	

Ø4.5mm




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Twist Drill (Ø3.0/4.1)	Cortical Drill (F4.5) Bottom line	Cortical Drill (F5.0) Top line	Ø4.5 Implant
Soft	▶		▶				Implant Placement
Normal	▶	▶		▶	▶		
Hard	▶	▶		▶		▶	

Ø5.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.1)	Twist Drill (Ø3.0/4.6)	Cortical Drill (F5.0) Bottom line	Cortical Drill (F5.0) Top line	Ø5.0 Implant
Soft	▶		▶				Implant Placement
Normal	▶	▶		▶	▶		
Hard	▶	▶		▶		▶	


Ø5.5mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Cortical Drill (F5.5) Bottom line	Cortical Drill (F5.5) Top line	Ø5.5 Implant
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	


Drilling Sequence **Ultra-wide 123 Twist Drill**  
**TSII Ultra-wide** | **SSII Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.2)	Twist Drill (Ø4.6/5.5)	Cortical Drill (F6.0)	Ø6.0 Implant
Soft	▶	▶	▶	▶			
Normal	▶	▶	▶		▶		Implant Placement
Hard	▶	▶	▶		▶	▶	


Ø7.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.5)	Twist Drill (Ø4.6/5.5 (F7.0 Soft))	Twist Drill (Ø5.5/6.5)	Cortical Drill (F7.0)	Ø7.0 Implant
Soft	▶	▶	▶	▶	▶			
Normal	▶	▶	▶	▶		▶		Implant Placement
Hard	▶	▶	▶	▶		▶	▶	


Drilling Sequence **Ultra-wide 123 Twist Drill**  
**TSIII Ultra-wide** | **SSIII Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.2)	Twist Drill (Ø4.6/5.5)	Cortical Drill (F6.0)	Ø6.0 Implant
Soft	▶	▶	▶	▶(6mm)			
Normal	▶	▶	▶		▶(6mm)		Implant Placement
Hard	▶	▶	▶		▶(6mm)	▶	

Ø7.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.5)	Twist Drill (Ø5.5/6.2 (F7.0 Soft))	Twist Drill (Ø5.5/6.5)	Cortical Drill (F7.0)	Ø7.0 Implant
Soft	▶	▶	▶	▶	▶(6mm)			
Normal	▶	▶	▶	▶		▶(6mm)		Implant Placement
Hard	▶	▶	▶	▶		▶(6mm)	▶	

Ultra KIT (HULTRK)


Applicable Products Ultra-wide

Top panel components

Open Wrench  
SPOW



Ratchet Wrench  
RCWC




OSSTEM<sup>®</sup>  
IMPLANT

• More details on KIT components can be found in Surgical Instruments(222p~242p)




Drilling Sequence **III Type** **Straight Drill**  
**TSIII | SSIII**  
(Length : 10mm)

Ø3.0mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Cortical Drill 2 (F3.0)	Ø3.0 Implant
Soft	▶			Implant Placement
Normal	▶	▶		
Hard	▶	▶	▶	

Ø3.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø2.2)	Cortical Drill 3 (F3.5)	Cortical Drill 3 (F3.5)	Ø3.5 Implant
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	


Ø4.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.3)	Cortical Drill 3 (F4.0)	Cortical Drill 3 (F4.0)	Ø4.0 Implant
Soft	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶		▶	


Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal/hard bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

Ø4.5mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Cortical Drill 3 (F4.5)	Cortical Drill 3 (F4.5)	Ø4.5 Implant
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

Ø5.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.3)	Cortical Drill 3 (F5.0)	Cortical Drill 3 (F5.0)	Ø5.0 Implant
Soft	▶	▶	▶	▶	▶				Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	▶		▶	


Ø5.5mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Cortical Drill 3 (F5.5)	Cortical Drill 3 (F5.5)	Ø5.5 Implant
Soft	▶	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	▶		▶	

Drilling Sequence **IV Type** **Straight Drill**  
**TSIV**  
(Length : 10mm)


Ø4.0mm



Ø4.5mm




Ø5.0mm




Drilling Sequence **Ultra-wide** **Straight Drill**  
**TSII Ultra-wide** | **SSII Ultra-wide**  
(Length : 10mm)

Ø6.0mm




Ø7.0mm




Recommended placement torque ≤ 40Ncm  
TSIV/USIV implants are specifically designed for sinus lift and soft bones, and not recommended for bone quality of normal bones or higher. In the case of TSIV/USIV implants, the speed of implant placement is fast due to the large pitch of the thread, and placing the implant by lowering the drilling speed to 15rpm or lower is recommended.

Drilling Sequence **Ultra-wide** **Straight Drill**  
**TSII Ultra-wide** | **SSII Ultra-wide**  
(Length : 10mm)


**Ø6.0mm**  
Drilling sequence with trephine in the healed mature bone



Immediate placement at the extraction socket




Immediate replacement of the failed implant




Drilling Sequence **Ultra-wide** **Straight Drill**  
**TSIII Ultra-wide** | **SSIII Ultra-wide**  
(Length : 10mm)

**Ø6.0mm**




**Ø7.0mm**




Recommended placement torque ≤ 40Ncm  
TS implant placement depth : For normal/hard bones, 1mm deeper than the bone level; for soft bones, matched to the bone level to maintain the stability

Drilling Sequence **Ultra-wide Straight Drill**  
**TSIV Ultra-wide**  
(Length : 10mm)

Ø6.0mm



Ø7.0mm





OSSTEM<sup>®</sup>  
IMPLANT

Recommended placement torque ≤ 40Ncm

485 KIT (O485K)

Applicable Products TSIII SSIII

Top panel components




OSSTEM<sup>®</sup>  
IMPLANT

• More details on KIT components can be found in Surgical Instruments(222p~242p)


Drilling Sequence **485 Drill**  
**TSIII | SSIII**  
(Length : 7mm)

Ø4.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.0)	485 Drill (F4.5)	Ø4.0 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	

Ø4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.5)	485 Drill (F5.0)	Ø4.5 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	

Ø5.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F5.0)	485 Drill (F5.5)	Ø 5.0 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	

OSSTEM<sup>®</sup>  
IMPLANT

Assist KIT (OAK) NEW 2025

- Bone profilers are only sold in the packing unit of “Guide Screw + Bone Profiler”
- For information on the order code for TS / KS Bone Profiler, please see page 402



OSSTEM<sup>®</sup>  
IMPLANT

• More details on KIT components can be found in Surgical Instruments(222p~242p)

Surgical Instruments

123 Guide Drill

- A drill for forming a hole to facilitate initial drilling
- Facilitating drilling depth adjustment by assembling a stopper


D	Ø2.0
OGD2027L	



Lance Drill (Guide Drill)

- Forming a hole to facilitate initial drilling
- Bone density determined through drilling


L	Short	Long
AGDSC		AGDLC



Sidecut Drill


- A drill to remove the side parts with the cutting edge of the body
- Used for removing the ridge of a fresh extraction socket
- Facilitating site preparation of a fresh extraction socket

L \ D1 / D2	Ø1.5 / 2.0	Ø2.0 / 2.5	Ø3.0 / 3.5
13	OSLMDS	OSLMD20S	-
16.5	-	-	OSLMD30L
20	OSLMDL	OSLMD20L	-



122 Taper Drill

- Included in 122 Taper KIT
- A dedicated taper drill for taper (III type) implants types available for each diameter and length
- Color coded handle indicates the implant diameter
- A drill slightly larger in diameter is used for removing cortical bone from hard bone
- F = Final drill




L \ TL	TL	F3.5	F4.0	F4.5	F5.0	F5.5	F6.0	F7.0
	Y-Dim.	0.7	0.9	1.0	1.0	1.0	1.0	1.0
4.0	29.5	122TPD3504	122TPD4004	122TPD4504	122TPD5004	122TPD5504	-	-
5.0	29.5	122TPD3505	122TPD4005	122TPD4505	122TPD5005	122TPD5505	-	-
6.0	30.5	122TPD3506	122TPD4006	122TPD4506	122TPD5006	122TPD5506	122TPD6006	122TPD7006
7.0	31.5	122TPD3507	122TPD4007	122TPD4507	122TPD5007	122TPD5507	122TPD6007	122TPD7007
8.5	33	122TPD3508	122TPD4008	122TPD4508	122TPD5008	122TPD5508	122TPD6008	122TPD7008
10	34.5	122TPD3510	122TPD4010	122TPD4510	122TPD5010	122TPD5510	122TPD6010	122TPD7010
11.5	34.5	122TPD3511	122TPD4011	122TPD4511	122TPD5011	122TPD5511	122TPD6011	122TPD7011
13	36	122TPD3513	122TPD4013	122TPD4513	122TPD5013	122TPD5513	122TPD6013	122TPD7013
15	38	122TPD3515	122TPD4015	122TPD4515	122TPD5015	122TPD5515	-	-
Color		Yellow	Green	Blue	Red	Yellow	Green	Blue

Surgical Instruments

Taper Drill


- Included in Taper KIT
- A dedicated taper drill for taper (III type) implants of each diameter and length
- A stopper drill with 1mm margin
- Color coded handle indicates the implant diameter
- F3.5: yellow, F4.0: green, F4.5: blue, F5.0: red, F5.5: yellow



L	TL	F3.5	F4.0	F4.5	F5.0	F5.5
	Y-Dim.	0.8	0.9	1.0	1.0	1.0
5.0	29.5	TPD3C3505	TPD3C4005	TPD3C4505	TPD3C5005	-
6.0	30.5	TPD3C3506	TPD3C4006	TPD3C4506	TPD3C5006	TPD3C5506
7.0	31.5	TPD3C3507	TPD3C4007	TPD3C4507	TPD3C5007	TPD3C5507
8.5	33	TPD3C3508	TPD3C4008	TPD3C4508	TPD3C5008	TPD3C5508
10	34.5	TPD3C3510	TPD3C4010	TPD3C4510	TPD3C5010	TPD3C5510
11.5	34.5	TPD3C3511	TPD3C4011	TPD3C4511	TPD3C5011	TPD3C5511
13	36	TPD3C3513	TPD3C4013	TPD3C4513	TPD3C5013	TPD3C5513
15	38	TPD3C3515	TPD3C4015	TPD3C4515	TPD3C5015	TPD3C5515
Color		Yellow	Green	Blue	Red	Yellow

Taper Ultra Drill


- Included in Taper Ultra KIT
- A dedicated taper drill for Taper Ultra-Wide implants of each diameter and length
- A stopper drill with 1mm margin
- Color coded handle indicates the implant diameter
- F = Final drill



L	F6.0	F7.0
6	TPD3C6006	TPD3C7006
7	TPD3C6007	TPD3C7007
8.5	TPD3C6008	TPD3C7008
10	TPD3C6010	TPD3C7010
11.5	TPD3C6011	TPD3C7011
13	TPD3C6013	TPD3C7013
Color		GreenBlue

485 Drill


- Included in 485 KIT
- A drill for placing short implants in alveolar bone lacking in vertical dimension
- Ø 2.2 drill : straight drill
- Except for Ø 2.2 drill, the top blade of the drill is in the shape of CAS Drill, and the side blade is in the shape of taper drill
- A stopper drill with 1mm margin
- Recommended drilling speed: 800~1,200rpm



L	Type	Ø2.2	Pilot	F4.0	F4.5	F5.0	F5.5
4.0		O485D2204	O485D3504	O485D4004	O485D4504	O485D5004	O485D5504
5.0		O485D2205	-	O485D4005	O485D4505	O485D5005	O485D5505
6.0		O485D2206	O485D3506	O485D4006	O485D4506	O485D5006	O485D5506
7.0		O485D2207	-	O485D4007	O485D4507	O485D5007	O485D5507
8.5		O485D2208	-	O485D4008	O485D4508	O485D5008	O485D5508

123 Twist Drill

- Included in 123 Straight Simple KIT
- A straight drill to reduce the number of drilling (marking drill)
- A color-coded handle of the 123 Drill indicates the drill diameter and the main implant used
- Facilitating drilling depth adjustment by assembling a stopper
- Use of a stopper is necessary because of the difficulty of controlling the depth due to excellent cutting force
- F = Final drill











TL	D1 / D2 / D3			
	F3.5(Ø2.2 / 3.0)	F4.0(Ø3.0 / 3.6)	F4.5(Ø3.0 / 3.6 / 4.1)	F5.0(Ø3.0 / 4.1 / 4.6)
34	2D2230FNS	2D3036FNS	2D3041FNS	2D3046FNS
40.4	2D2230FNL	2D3036FNL	2D3041FNL	2D3046FNL
Color		Yellow	Green	Blue

Surgical Instruments


123 Drill Stopper

- Included in 123 Straight Simple KIT
- Number on the stopper indicates the protruding length of the tip when assembled to a drill or other instruments
- Color coded by length to facilitate estimation of the length and repositioning of the KIT

L	6.2	7	8	9.5	11	12.5	14	16
								
	ODST05	ODST06	ODST07	ODST08	ODST10	ODST11	ODST13	ODST15

123 Ultra Twist Drill


- Included in 123 Straight KIT/123 Straight Full KIT
- A 2-stage drill with functions of both Pilot Drill and Twist Drill
- A straight drill to reduce the number of drilling (with stopper)
- A dedicated drill is used for F7.0 implant in soft bone
- F = Final drill



L	TL	F3.5(Ø4.6 / 5.2)	F6.0(Ø4.6 / 5.5)	F7.0Soft(Ø5.5 / 6.2)	F7.0(Ø5.5 / 6.5)
6	30.5	3D465206T	3D465506T	-	3D556506T
7	31.5	3D465207T	3D465507T	-	3D556507T
8.5	33.5	3D465208T	3D465508T	-	3D556508T
10	34.5	3D465210T	3D465510T	-	3D556510T
11.5	34.5	3D465211T	3D465511T	-	3D556511T
13	36.0	3D465213T	3D465513T	3D556213T	3D556513T
Color		Yellow	Green	Blue	Red

123 Twist Drill (Stopper Drill)


- Included in 123 Straight KIT/ 123 Straight Full KIT
- A straight drill to reduce the number of drilling (with stopper)
- A color-coded handle of the 123 Drill indicates the drill diameter and the main implant used
- F = Final drill



		D1 / D2 / D3			
L	TL	F3.5(Ø2.2 / 3.0)	F4.0(Ø3.0 / 3.6)	F4.5(Ø3.0 / 3.6 / 4.1)	F5.0(Ø3.0 / 4.1 / 4.6)
6	30.5	2D223006LC	2D303606LC	2D304106LC	2D304606LC
7	31.5	2D223007LC	2D303607LC	2D304107LC	2D304607LC
8.5	33	2D223008LC	2D303608LC	2D304108LC	2D304608LC
10	34.5	2D223010LC	2D303610LC	2D304110LC	2D304610LC
11.5	34.5	2D223011LC	2D303611LC	2D304111LC	2D304611LC
13	36	2D223013LC	2D303613LC	2D304113LC	2D304613LC
15	38	2D223015LC	2D303615LC	2D304115LC	2D304615LC
Color		Yellow	Green	Blue	Red

Twist Drill (Stopper Drill)

- Included in New Hanaro KIT
- Long stopper (6mm)  
: Enabling a procedure without drill extension for posterior region
- The color coded stopper indicates the drill length.




L	TL	D	Ø2.2	Ø3.0	Ø3.3	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
	Y-Dim		0.6	0.9	1.0	1.0	1.0	1.0	1.0	1.0
6	30.5	2D2206LC	3D3006LC	-	-	3D3806LC	-	-	-	-
7	31.5	2D2207LC01	3D3007LC01	-	-	3D3807LC01	-	-	-	-
8.5	33	2D2208LC01	3D3008LC01	-	-	3D3808LC01	-	-	-	-
10	34.5	2D2210LC01	3D3010LC01	-	-	3D3810LC01	-	-	-	-
11.5	34.5	2D2211LC01	3D3011LC01	3D3311LC01	3D3611LC01	3D3811LC01	3D4111LC01	3D4311LC01	3D4611LC01	-
13	36	2D2213LC01	3D3013LC01	-	-	3D3813LC01	-	-	-	-

Surgical Instruments

Twist Drill (Non Stopper Drill)

- Included in New Hanaro KIT
- Used for limited access for the Stopper Drill into the oral cavity
- See the image provided in the Non-stopper Drill section for the sizes of the drill marking lines for short/long types




TL \ D	Ø1.5	Ø2.0	Ø2.2	Ø2.7	Ø3.0	Ø3.3
33	2D1518FNLC	2D2018FNLC	2D2218FNLC	3D2718FNLC	3D3018FNLC	3D3318FNLC
41	-	-	2D2215FNLC01	3D2715FNLC01	3D3015FNLC01	3D3315FNLC01

TL \ D	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
33	3D3618FNLC	3D3818FNLC	3D4118FNLC	3D4318FNLC	3D4618FNLC
41	3D3615FNLC01	3D3815FNLC01	3D4115FNLC01	3D4315FNLC01	3D4615FNLC01

Direct Drill

- Included in Ultra KIT
- A 2-stage drill with functions of both Pilot Drill and Twist Drill
- Enabling final drilling without pilot drilling
- Increased primary stability in a fresh extraction socket with reduced dead space in the apex




D1 / D2	Ø4.6 / 5.2	Ø4.6 / 5.5	Ø5.5 / 6.2	Ø5.5 / 6.5
	3D5213FNLC	3D5513FNLC	3D6213FNLC	3D6513FNLC

Long Shank Pilot Drill

- Included in New Hanaro KIT
- Used for adjusting the path of the drilling hole
- Maintains the previous drilling path for the next drilling


D1 / D2	Ø2.0 / 2.7	Ø2.0 / 3.0	Ø3.0 / 3.8	Ø3.0 / 4.1
	APD270C	APD300C	APD380C	APD410C



Taper Cortical Drill  
(Taper Implant TSIII, SSIII)

- Included in Taper KIT
- A drill used for removing cortical bone from hard bone (used right after the use of Taper Drill)
- A dedicated drill equipped for each implant diameter
- F3.5~5.0 drill marking line : bottom line for placing implants of 8.5mm or smaller, and top line for implants of 10mm or larger
- F5.5 drill marking line : bottom line for placing implants of 6mm or smaller, midline for 7mm implants, and top line for implants of 8.5mm or larger
- Drilling up to the lower marking line recommended
- F = Final drill


F3.5	F4.0	F4.5	F5.0	F5.5
TCD4C35	TCD4C40	TCD4C45	TCD4C50	TCD4C55



Cortical Drill (Ultra-Wide)

- Included in Ultra KIT, 122 Taper / Taper KIT (for ultra-wide)
- A drill used for removing cortical bone from hard bone (for ultra-wide)
- A dedicated drill equipped for each implant diameter
- Drilling up to the lower marking line recommended
- F = Final drill


F6.0	F7.0
CD4C60	CD4C70



Surgical Instruments

123 Cortical Drill

- Included in 123 Straight Simple KIT, 123 Straight KIT and 123 Straight Full KIT
- A drill used for removing cortical bone from hard bone
- Drilling up to the lower marking line recommended
- II type marking line : for hard bone
- III type marking line : bottom line for normal bone, and top line for hard bone
- IV type marking line : for normal bone
- A color-coded handle of the 123 Drill indicates the drill diameter and the main implant used
- F = Final drill




Type	F3.5	F4.0	F4.5	F5.0
II	O2CD35	O2CD40	O2CD45	O2CD50
III	O3CD35	O3CD40	O3CD45	O3CD50
Color	Yellow	Green	Blue	Red

Tapered Implant Tap  
(Taper Implant TSIII, SSIII)

- A dedicated tap for Tapered Implant (III type)
- Used for hard bones, forming the implant thread shape
- A torque wrench is used after connecting to the engine (25rpm recommended) or a mount extension
- Tapping up to the bottom marking line recommended (For F5.0, the bottom line for placing 7.0mm or smaller implants, and top line for 8.5mm or greater implants)
- F = Final drill

F3.5	F4.0	F4.5	F5.0
OFTS35	OFTS40	OFTS45	OFTS50



Cortical Drill 2 (TSII, SSII)

- Included in New Hanaro KIT
- A drill used for removing cortical bone from hard bone (For II type)
- A dedicated drill equipped for each implant diameter
- Drilling up to the lower marking line recommended
- F = Final drill


F3.5	F4.0	F4.5	F5.0
CD2C35	CD2C40	CD2C45	CD2C50



Parallel Pin (122 Taper Drill)

- Included in 122 Taper KIT
- A dedicated parallel pin for 122 Taper Drill
- Used for checking the position and direction of bone preparation
- Bottom part for the Ø2.2 drill, and top part for the guide drill

APP2227
---------



Cortical Drill 3  
(Taper Implant TSIII, SSIII)

- Included in New Hanaro KIT
- A drill used for expanding the cortical bone after the use of Straight Drill
- Used after forming the final drill hole in normal or harder bone
- A dedicated drill equipped for each implant diameter
- Marking line: bottom line for normal bone, and top line for hard bone
- Drilling up to the lower marking line recommended


F3.0	F3.5	F4.0	F4.5	F5.0	F5.5
CD4C30	CD4C35	CD4C40	CD4C45	CD4C50	CD4C55



Parallel Pin (Taper Drill)

- Included in 122 Taper KIT and Taper KIT
- A dedicated parallel pin for Taper Drill
- Used for checking the position and direction of bone preparation
- Bottom part for implant diameter drill, top part for Initial Drill
- Color coded according to the implant diameter (F3.5 : yellow, F4.0 : green, F4.5 : blue, F5.0 : silver)

F3.5	F4.0	F4.5	F5.0
TPP3522	TPP4022	TPP4522	TPP5022




Surgical Instruments

Parallel Pin (123 Drill)


- Included in 123 Straight Simple KIT, 123 Straight KIT and 123 Straight Full KIT
- Used for checking the position and direction of bone preparation
- Bottom part for Initial Drill, top part for F3.5(Ø2.2/3.0) Drill

D	Ø4.0	Ø5.0
	OPLP400	OPLP500



Tissue Punch

- Included in Assist KIT
- Instrument used for flapless surgery
- Marked at 2mm intervals for measuring gingival height
- Packing unit : tissue punch + guide pin
- ※ Using Tissue Punch with a diameter smaller by 0.7~1.5mm than the Healing Abutment is recommended
- Recommended drilling speed: 1,000~ 1,200rpm




D	Ø3.3	Ø3.8	Ø4.3	Ø4.8	Ø5.3
	OSTP33	OSTP38	OSTP43	OSTP48	OSTP53
TS	Ø 4.0/4.5	Ø 4.5/5.0	Ø 5.0	Ø 6.0	Ø 6.0
SS	-	Ø 4.8	-	Ø 6.0	Ø 6.0
US	Ø 4.0	Ø 5.0	Ø 5.0	Ø 6.0	Ø 6.0

For healing abutment applications

Parallel Pin

- Included in New Hanaro KIT
- Used for checking the position and direction of bone preparation


D	Ø4.0	Ø5.0	Ø6.0	Full Set
	APP400	APP500	APP600	APPS



Trial Pin (Ultra-wide)

- 123 straight full kit, Included in Ultra KIT
- Checking the internal width and depth of a fresh extraction socket or a failed implant socket
- Checking the drilling depth after using Direct Drill as the final drill
- Used as Parallel Pin

D	Ø5.2	Ø5.5	Ø6.2	Ø6.5
	UWFTP52	UWFTP55	UWFTP62	UWFTP65



Bone Profiler (TS)

- Included in Assist KIT
- Used for removing bone around the implant for the 1st and 2nd stage surgery
- Used by connecting Guide Screw to the implant and removing the bone to adjust the shape of the Healing Abutment
- Guide Screw protecting the morse taper entrance of the implant
- Packing unit : bone profiler + guide screw
- Recommended drilling speed: 50rpm
- C = Connection




C \ D (Healing Abutment)	Ø4.0	Ø4.5	Ø5.0	Ø6.0	Ø7.0
TS Mini/Regular	GSBP40 Mini+Regular guide screw	GSBP45 Mini+Regular guide screw	GSBP50 Regular guide screw	GSBP60 Regular guide screw	GSBP70 Regular guide screw

Surgical Instruments

Trephine Drill


- Included in Assist KIT
- Used for bone collection or for removing damaged or failed implants
- Used for removing septal bone
- Used as an Initial Drill for Ultra-wide implants
- Recommended drilling speed: 1,000~1,200rpm



L \ D (Inner / Outer)	3.7 / 4.5	4.2 / 5.0	4.7 / 5.5	5.2 / 6.0	5.7 / 6.5	6.2 / 7.0
Short	TD37S	TD42S	TD47S	TD52S	TD57S	TD62S
Long	TD37	TD42	TD47	TD52	TD57	TD62

Crest Remover


- Included in Assist KIT
- Marking the implant placement position after removing the narrow alveolar ridge horizontally
- Recommended drilling speed
  - Angled type : 1,200~1,500rpm
  - Straight type : 15,000~30,000rpm



L \ D	Ø5.0	Ø7.0
29	CERM50A	CERM70A
45	CERM50S	-

AutoBone Collector®

- Included in Assist KIT
  - Used for autogenous bone collecting
  - Comes in a Drill + Stopper set
  - Recommended drilling speed : 300~600rpm
  - Number of uses for the drill and stopper: 50 times
- ※ Before initial drilling, connect the stopper to the first stage locking and harvest autogenous bone while drilling 4mm into the second stage locking (after harvesting, stop the drill and remove it as is with autogenous bone kept in the stopper)



L \ D	Ø3.0	Ø4.0	Ø5.0	Ø6.0
Short (18.94)	ABC304S	ABC404S	ABC504S	ABC604S
Long (21.94)	ABC304L	ABC404L	ABC504L	ABC604L

Drill Extension

- Extending the length of a drill or other hand piece tool (drill extended by 16.9mm)
- Risk of bending or fracture upon exerting excessive force on inadequate assembly
- Common component of Taper KIT and Straight KIT


L (Extention)	16.9
	ODE



TS NoMount Driver

- Driver directly connected to the fixture upon placing with a surgical hand piece
- C = Connection


L \ C	Mini	Regular
Short	TSNMDMS	TSNMDRS
Long	TSNMDML	TSNMDRL
Ex.Long	TSNMDME	TSNMDRE



SS NoMount Driver

- Driver directly connected to the fixture upon placing with a surgical hand piece
- C = Connection

L \ C	Regular / Wide
Short	SSNMDS
Long	SSNMDL
Ex.Long	SSNMDE




Surgical Instruments

TS NoMount Torque Driver

- Driver directly connected to the implant upon placing with a wrench
- Make sure to check and confirm proper assembly before use  
(Risk of fracture even at low torque in case of inadequate assembly)
- Note that it cannot be removed in case of fracture
- C = Connection


L \ C	Mini	Regular
Short	GSNMT32S	GSNMT35S
Long	GSNMT32L	GSNMT35L
Ex.Long	GSNMT32E	GSNMT35E



SS Implant Driver

- Used by assembling directly to the implant for final placement depth adjustment or removal
- C= Connection


L \ C	Regular / Wide
Short	SSRFDS
Long	SSRFDL
Ex.Long	SSRFDE



SS NoMount Torque Driver

- Driver directly connected to the implant upon placing with a wrench
- Make sure to check and confirm proper assembly before use  
(Risk of fracture even at low torque in case of inadequate assembly)
- Note that it cannot be removed in case of fracture
- C = Connection


L \ C	Regular / Wide
Short	SSNMT39S
Long	SSNMT39L



Simple Mount Driver

- Used by assembling to the simple mount for implant placement


L	
Short	ASMDS
Long	ASMDL



Simple Mount Extension

- Used by connecting to a wrench for extending the simple mount length or applying torque manually


L	
Short	ASMES
Long	ASMEL



TS Implant Driver

- Used by assembling directly to the implant for final placement depth adjustment or removal
- C = Connection

L \ C	Mini	Regular
Short	GSMFDS	GSRFDS
Long	GSMFDL	GSRFDL
Ex.Long	GSMFDE	GSRFDE



Torque Extension

- Extending the length of the instrument used by connecting to a wrench (extension by 10mm)

	OTE
--	-----




# Surgical Instruments

## Removal Tool (Implant Mount)

- Used after removing mount screw in case of jamming between the implant and mount
- Used by assembling to the driver handle and torque wrench
- Removing mount by rotating forwarding after inserting vertically
- App = Application


App	Mini (TS)	Regular (TS,SS)	Wide (SS)
	ERFM	HRFR	



## Positioning Guide

- Instrument facilitating drilling interval setting for implant placement
- Placed into the hole for use after initial drilling
- Packing unit : individual component or in a set

W/L	2.5 / 21.5	6.0 / 17.5	11 / 17.5
	APG201	APG202	APG203



## Tissue Height Gauge (TS)

- Instrument to measure the gingival height by assembling to the implant connection for top G/H selection in TS implant placement

HGTSHG
--------



## Depth Gauge

- Used for measuring the drilling depth (7-15mm) or as an open wrench
- A common component of 122 Taper & Taper KIT

OSDG
------



## Simple Open Wrench

- Used for removing a simple mount from weak bone
- Easy intraoral placement with a neck angle of 30°

ASOW
------



## Ratchet Wrench

- A dedicated wrench for anti-reverse procedure
- Excessive torque exertion may result in internal damage to bone or implant


CITQW-1185A
-------------



## L-Wrench

- 1.2 hex driver for overcoming narrow spacing
- Torque indication : 5~8Ncm torque at the level when the wrench appears to be bent a little (within 10°)

LWC
-----



## Torque Wrench (Spring Type)

- Wrench to apply a constant torque (10/20/30Ncm) to screws and abutments
- When the set torque is applied, the neck of the Torque Wrench is bent for indication
- If a continuous force is applied while the neck is bent, it will cause application of excessive torque, resulting in screw fracture.


TW30
------



# Surgical Instruments

## Torque Wrench (Bar Type)


- Used for adjusting the implant placement position and tightening screws and abutments
- Applying torque according to the line marked with the torque value to be applied by pulling the bar



TW30B

## Torque Wrench Set


- A set of a two-way Torque Wrench and a Torque Connector
  - Applying forward/reverse torque by rotating the Torque Wrench handle without removing the connector
- Compatible with the Machine Driver connector of OSSTEM
- Applying torque according to the line marked with the torque value to be applied by pulling the bar
- Packing unit : changeable torque wrench + torque connector



MX30

## Bone Mill

- Forming particulate bone with collected autogenous bone



ABM

## Anterior Hand Driver (Implant)

- Instrument for manual placement in anterior region
- Used by connecting to a NoMount Torque Driver or Implant Driver
- Excessive torque may result in fracture of the implant or driver



AHDI

## Torque Handle


- Manual instrument used by connecting to the contra-angle hand piece (1:1 gear ratio for hand piece)
- Used for tightening screws such as Healing Abutment, Cover Screw, Abutment Screw, and Orthodontic Screw (used for temporary tightening of Abutment Screw, which requires final tightening with a Torque Wrench)
- Excessive torque may result in fracture or malfunction of the hand piece



TQHD

## Torque Connector


- Connector for connecting the torque square driver with a two-way Torque Wrench



ORC

## Machine Driver Connector


- Connector for connecting the Machine Driver with a two-way Torque Wrench



OMDC

## Driver Handle

- Used by connecting to the Torque Driver



TIDHC

# Surgical Instruments

## Machine Driver Handle


- Enabling hand rotation by connecting to any surgical instrument for engine



OMDH

## Finishing Reamer Set


- Used for removing lip from the inside of the casted body after casting plastic coping



FRSC

### Reamer user guide

1. Connected to the casted burn-out cylinder by selecting the reamer tip of the same size as the abutment
2. Rotating the reamer bite with constant force by holding the casted body
3. Reaming until no cutting occurs



OSSTEM<sup>®</sup>  
IMPLANT

Prosthetic Simple KIT (OPSK)

Prosthetic KIT (OPK)


Top panel components

Torque Wrench  
TW30B




SURGICAL KIT

244




SURGICAL KIT

245



Hand Driver


- Manual driver
- Tip holding feature (except internal hex type)
- Internal hex type length:11



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Ex.Short (8)	-	AHD09MSH	AHD12MSH	-	-
Short (13)	ASD05SH	AHD09SH	AHD12SH	IHD20H	IHD27H
Middle (15)	-	-	AHD12MH	-	-
Long (18)	ASD05LH	AHD09LH	AHD12LH	-	-
Ex.Long (25)	-	-	AHD12EH	-	-

Torque Driver


- A driver for Torque Wrench assembly
- Tip holding feature
- Compliance to the recommended torque is required(excessive torque may result in fracture)
- Risk of fracture even at low torque in case of inadequate assembly
- Exerting torque with the driver straight up is required (with no tilting)
- Be sure to replace any bent tips due to extended use or excessive torque application



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Ex.Short (8)	-	-	TRHD12MS	-	-
Short (13)	TRSD05S	TRHD09S	TRHD12S	TIHD20S	-
Middle (15)	-	-	TRHD12M	-	-
Long (20)	TRSD05L	TRHD09L	TRHD12L	TIHD20L	TIHD27
Ex.Long (25)	TRSD05E	-	TRHD12E	-	-

Machine Screw Driver

- Driver for engine
- Tip holding feature (except internal hex type)
- Internal hex type length : 8



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Osstem Torque(5)	-	-	OTH12S	-	-
Short (5.6)	AMSD05S	AMSD09S	AMSD12S	-	-
Long (11.6)	AMSD05L	AMSD09L	AMSD12L	EIHD20	EIHD27
Ex.Long (17.6)	-	-	AMSD12E	-	-


Application

For driver applications  
(Commonly used for Hand, Machine Screw, and Torque Drivers)

Cover screw (US mini)	Healing abutment, Gold/NP-Cast abutment, Cemented abutment screw, Mount screw	Esthetic abutment screw Regular, Esthetic-low abutment screw, standard	Wide esthetic-low abutment screw
-----------------------	---	--	----------------------------------

Angled Torque Driver

- A driver for Torque Wrench assembly
- No holding feature
- Recommended tightening torque: 30Ncm (excessive torque may result in fracture)
- Do not remove the tube preventing debris upon fracture
- Recommended number of use cycles: 10 times
- Set : 3ea




L \ Type	1.2 Hex	1.2 Hex (Set)
Short (13)	ATD12S	ATD12S3S
Long (20)	ATD12L	ATD12L3S

Repair Torque Driver


- Handle diameter reduced compared to Torque Driver (Ø2.1 → Ø1.6)
- Enables minimization of crown hole diameter for prosthesis repair or SCRP procedure

L \ Type	1.2 Hex
Short (13)	TRHD12SR
Long (20)	TRHD12LR



Solid Abutment Driver

- A dedicated driver for Solid Abutment
- Applying torque after inserting the groove of the Solid Abutment to the part of the driver with a triangular marking
- Recommended tightening torque: 30Ncm



Regular

L \ Type	Square	Round
Short (6)	SDSS	SDRS
Long (12)	SDSL	SDRL

Wide

L \ Type	Square
Short (10)	SD60S

Excellent Solid Abutment Driver

- A dedicated driver for Excellent Solid Abutment
- Applying torque after inserting the groove of the Excellent Solid Abutment to the part of the driver with a triangular marking
- Recommended tightening torque: 30Ncm



Regular

L \ Type	Square	Round
Short (6)	ESDSS	ESDRS
Long (12)	ESDSL	ESDRL


Wide

L \ Type	Square
Short (10)	ESD60S

O-ring Abutment Driver

- A dedicated driver for O-ring Abutment

	Small	Normal
	STAOD	AORD



Rigid Outer Driver

- A dedicated driver for Rigid Abutment
- Recommended tightening torque: 30Ncm


L \ D (Abutment)	Ø4.0	Ø4.5	Ø5.0	Ø6.0
Short (16.5)	ORDMS	ORD45S	ORDRS	ORDWS
Long (21.5)	ORDML	ORD45L	ORDRL	ORDWL



Octa Abutment Driver

- A dedicated driver for Octa Abutment
- Recommended tightening torque: 30Ncm

L \ Type	Square	Round
Short	ODSS	ODRS
Long	ODSL	ODRL



Multi Abutment Machine Driver

- A dedicated Machine Driver for Multi Abutment



Abutment Holder

- Supplementary instrument for convenient connection of a 2-piece abutment which is difficult to hold with a hand in all intraoral regions



Abutment Positioning Driver

- Used for assembling the abutment in the prosthetic process after placing an implant
  - ※ For Transfer Abutment only
- Function to help convenient and stable mounting and tightening of the abutment being pushed away by gingiva
- Used according to the H and G/H lengths of the abutment to be removed as shown below

(Unit : mm)

Range of use	Short					Long				
	≤9					≥10				
H + G/H	5	6	7	8	9	10	11	12	13	14



Multi Abutment Outer Driver

- A dedicated Torque Driver for Multi Abutment




Locator® Torque Driver

- A dedicated Torque Driver for Locator Abutment



Osstem Torque Driver


- Dedicated driver for Osstem Torque, which may not be compatible (connected or separated) with a general hand piece
- Used after matching the triangle on the outside of the driver with the groove or side of the abutment
- Solid and Excellent Solid Driver are only compatible with Ø4.8
- 1.2 hex type L : 5



L \ Type	1.2 Hex	Rigid 4.0	Rigid 4.5	Rigid 5.0	Rigid 6.0	Solid	Excellent Solid
Short(10)	OTH12S	OTR40S	OTR45S	OTR50S	OTR60S	OTS48S	OTE48S
NEW 2025 Long(15)	-	OTR40L	OTR45L	OTR50L	OTR60L	OTS48L	OTE48L

Path Probe (TS)

- Checking the path and measuring the gingival height after TS Implant placement
- C = Connection



Reamer Bite

- Cutting edge to remove lip from the inside of the casted body after casting plastic coping




FRBC

Reamer Tip (Rigid Abutment)

- Guide part inserted into the casted body for removing lip from the inside after casting plastic coping (for Rigid Abutment)


D	Ø4.0	Ø4.5	Ø5.0	Ø6.0
	GSRFRT400	GSRFRT450	GSRFRT500	GSRFRT600



Reamer Tip (Solid, Excellent Solid Abutment)

- Guide part inserted into the casted body for removing lip from the inside after casting plastic coping
- For Solid Ø6.0 and Excellent Solid Ø4.8
- P= Platform

P	Regular(Ø4.8)	Wide(Ø6.0)
Solid	FRTS480	FRTS600
Ex.Solid	FRTE480	FRTE600




CAS KIT (HCRSNK)

CAS KIT Surgical Instruments

Applicable Products TSIII / IV SSIII

Top panel components Lower panel components



Stopper (7)

SNST7 Yellow

Stopper (2,8)

SNST2 SNST8 Purple

Guide Drill

SNGD2027TL

Twist Drill

SNTD2213TS SNTD2213TL

Hydraulic Membrane Lifter

OCHML

Stopper (6,12)

SNST6 SNST12 Blue

Stopper (5,11)

SNST5 SNST11 Purple

Stopper (3,9)


SNST3 SNST9 Blue

Stopper (4,10)

SNST4 SNST10 Yellow

CAS Drill

- Safe lifting of the membrane while forming conical bone for maxillary sinus lift procedure
- Excellent bone removal at low-high speed, and collection of autogenous bone at low speed
- Stopper assembled for safe lifting
- Diameter of Final Drill is selected based on the bone quality regardless of Straight or Tapered Implant type
- Recommended drilling speed: 400~800rpm (400rpm for first use)




L \ D	Ø2.8	Ø3.1	Ø3.3	Ø3.6	Ø3.8	Ø4.1
Short	SNDR2813TS	SNDR3113TS	SNDR3313TS	SNDR3613TS	SNDR3813TS	SNDR4113TS
Long	SNDR2813TL	SNDR3113TL	SNDR3313TL	SNDR3613TL	SNDR3813TL	SNDR4113TL

Guide Drill

- Drill for marking of the implant placement position
- Used for removing side walls in a fresh extraction socket with formation of side edges
- Marking line at 2mm from the tip


L \ D	Ø2.0 / 2.7
Short	SNGD2027TL
Long	



Twist Drill (Ø2.2)

- Drilling 1mm under the remaining bone recommended
- Stopper assembled for safe lifting
- End line tip: 0.6mm


L \ D	Ø2.2
Short	SNTD2213TS
Long	SNTD2213TL



• For ordering codes of single item of CAS KIT, see 255~258 Page


Hydraulic Membrane Lifter Set

- Instrument for hydraulic lifting of sinus membrane
- Winged design with optimized sealing



Bone Condenser

- Instrument to push in the bone material into the sinus




Stopper

- Number on the stopper indicates the protruding length of the tip when assembled to a drill or other instruments
- Color coded by length
- Number of uses for the drill and stopper: 50 times



Hydraulic Membrane Lifter Tube

- Connected to the hydraulic membrane lifter




Bone Carrier

- Used for filling the inside of the sinus with bone
- Mounting the head by fastening the back of the body
- Replaceable head (SNBCH30 or SNBCH35) for use




Bone Carrier Head

- Used for filling the inside of the sinus with bone
- SNBCH30 : used after drilling with CAS Drill Ø3.1/3.3
- SNBCH35 : used after drilling with CAS Drill Ø3.6/3.8/4.1 drilling
- Used repeatedly by filling the back of the marking line of the head and taking little by little with a bone condenser to completely fill the inside of the sinus



Membrane Lifter

- Safe lifting of the membrane due to the round shape with no cutting edge
- Lifter selected according to the CAS-Drill diameter as membrane lifting is performed after using the CAS-Drill (head diameter is CAS Drill diameter - 0.2mm)
- CAS Stopper assembled and used for adjusting the depth
- Recommended drilling speed: 400~800rpm (400rpm for first use)
- Make sure to use a drill with irrigation




L \ D	Ø2.6	Ø2.9	Ø3.1	Ø3.4	Ø3.6	Ø3.9
Short	SNML2813TS	SNML3113TS	SNML3313TS	SNML3613TS	SNML3813TS	SNML4113TS
Long	SNML2813TL	SNML3113TL	SNML3313TL	SNML3613TL	SNML3813TL	SNML4113TL

Depth Gauge

- For checking internal lifting of the sinus and measuring the remaining bone depth


SNDG



Bone Spreader

- Instrument for spreading the filled bone using the engine
- Assembled with a stopper for use
- Recommended drilling speed: ≤30rpm (low speed mode)

D	Ø2.0	Ø3.0
	SNBS2015T	SNBS3015T



Y-Connector


- Y-shaped connector for hydraulic lifting of two drilling holes at the same time

SNYCT



LAS KIT (HLRSNK)


- Lateral Approach - Sinus KIT (LAS KIT) : KIT optimized for lateral approach in maxillary sinus lift procedure
- Including dome drill and core drill for safe formation of a lateral window; and Ø5.5/7.0 diameters available according to the size of the window
- Depth can be adjusted by mounting a stopper on the LAS Drill, and the window can be safely formed without perforating the membrane




• For ordering codes of single items of LAS KIT, see pages 262-263

LAS Full KIT (HLRSNKP)

- KIT with 6 additional sinus lift instruments to LAS KIT



LAS KIT Plus Lower Plate




• For ordering codes of single items of LAS KIT, see pages 262-263  
• Lower panel components of LAS Full KIT are the same as those of Sinus KIT (See page 310)

Dome Drill

- Forming a window while collecting bone
- Enhanced cutting force with macro and micro cutting edges in combination
- Depth adjusted by assembling with a stopper
- Recommended drilling speed: 1,200~1,500rpm
- ※ Over-drilling may result in damage to the membrane

L \ D	Ø5.5	Ø7.0	Wide Ø7.0
25	LSDR554TD	LSDR74TD	LSDR74WTD



Core Drill


- Forming a window while forming the bone lid
- Excellent cutting force and membrane stability owing to the CAS Drill design concept
- Recommended drilling speed: 1,200~1,500rpm
- ※ Over-drilling may result in damage to the membrane






L \ D	Ø5.5	Ø7.0
25	LSDR554TC	LSDR74TC




Side Wall Drill

- Expanding the window after drilling with a dome drill
- Cutting at 1mm above the bottom of the drill edge recommended
- Recommended drilling speed: 1,500rpm




Side Cutting edge height (mm)	1.0	2.0	3.0	4.0	5.0
CAS KIT stopper (mm)	8.0	9.0	10	11	12
Side wall drill + CAS KIT stopper					



※ Depth adjusted by the use of CAS KIT Stopper in common

Bone Separator

- Removing the bone lid from the inside of the core drill



Stopper


- Number on the stopper indicates the protruding length of the tip when assembled to a drill or other instruments
- Color coded by length
- Number of uses for the drill and stopper: 50 times

L	0.5	1.0	1.5	2.0	2.5	3.0
						
	LSNST0.5	LSNST1.0	LSNST1.5	LSNST2.0	LSNST2.5	LSNST3.0
Color	Yellow	Purple	Blue	Yellow	Purple	Blue

Denture 4U KIT (OD4UK)

Denture 4U KIT Surgical Instruments

Applicable Products **TSII / III**



Denture 4U Guide


- Guide for stable and accurate initial and intermediate drilling for Denture 4U procedure
- Anterior guide : drilling positioning for Ø2.2 in anterior region (tooth number 2 and 3 positions marked)
- Posterior guide : drilling positioning function for Ø3.0 drill in posterior region
- ※ Used by assembling with the posterior guide of desired angle
- A removable handle for Denture 4U Guide



D4UG

Posterior Guide

- Used by assembling to the anterior guide prior to procedure
- ※ Assembled with the angle marking side shown
- Enables adjustment of the implant placement position in posterior region and buccolingual inclination
- Prior to procedure, selecting the angle of Posterior Guide through CT scans is recommended
- ※ Replaceable during procedure
- Drilling is performed by slowly entering the guide hole, referring to the marking line on the side of the posterior guide hole
- Drilling depth adjusted by drilling to the bottom of the marking line in the mesial direction
- Marking line spacing on the rod : 2mm




Bottom of the marking line → Check mesial direction

Degree	0°	17°	30°
	D4UPG0	D4UPG17	D4UPG30

Crest Remover

- Used for bone flattening for Denture 4U Guide procedure
- Marking the implant placement position after removing narrowed ridge horizontally
- Recommended drilling speed
  - Angled type : 1,200~1,500rpm
  - Straight type : 15,000~30,000rpm



L	D	Ø5.0
29		CERM50A
45		CERM50S

Denture 4U KIT Surgical Instruments

Anchor Screw

- Used for stable fixing of the bone in place by connection to the fixed center hole of the Denture 4U Guide and the fixed hole of Posterior Guide
- Fixing the Anchor Screw with the Mount Driver; if the Anchor Screw is not fixed well at this time, the Anchor Drill should be used first for drilling
  - ※ Anchor drill used first for normal/hard bone
- Selecting an Anchor Screw of appropriate length according to the degree of posterior segment retraction
- Engine is stopped to prevent Anchor Screw from spinning with no traction when in contact with the guide

L \ D	Ø1.65
5	D4UAS5
8	D4UAS8
11	D4UAS11



Anchor Drill

- Used to form a hole in normal/hard bone prior to tightening an Anchor Screw
- Drilling with a 3mm drill is recommended prior to additional drilling with an 11mm drill

L \ D	Ø1.65
3	D4UAD3
11	D4UAD11



Guided Initial Drill

- Used for drilling in anterior region : Ø2.2 drilling into the anterior guide hole of the Denture 4U Guide
- Drilling is performed by selecting a desired drilling hole of the Anterior Guide
- Recommended drilling speed: 800rpm

L \ D	Ø2.2
5	GD2208NC



Twist Drill


- Drilling is performed by slowly entering the guide hole, with the angle matched as much as possible, referring to the marking line on the side of the Posterior Guide hole
- Drilling depth adjusted by drilling to the bottom marking line in the mesial direction
- Marking line spacing of the rod: 2mm
- Recommended drilling speed: 800rpm

L \ D	Ø3.0
5	D4U2D3005
18	D4U2D3018



Conutersink

- Drill for using the Taper Drill after removing the Denture 4U Guide
  - ※ For removing bone interference from the stopper of the Taper Drill
- Removing bone interference upon assembling to the Multi Angled Abutment




Indicator

- Checking the location of the mental foramen, and the placement direction and length of the implant beforehand to ensure the stability in the procedure
  - ※ For checking the location of the mental foramen by opening a flap completely



Path Checker

- Checking the location of the mental foramen by predicting the extended line of the path checker through panoramic or CT scan
  - ※ For checking the location of the mental foramen without fully opening a flap



# Denture 4U KIT

Surgical Instruments

## Simple Mount Driver

- Used for placing an Anchor Screw for stable fixing of the Denture 4U Guide in place

L  
Short ASMDS



## Multi Abutment Machine Driver

- A dedicated Machine Driver for Multi Abutment

MAMD



## Multi Abutment Outer Driver

- A dedicated Torque Driver for Multi Abutment

MAOD




Positioning Guide KIT (OPGPK)

Positioning Guide Full KIT (OPGAK)


SURGICAL KIT

270



SURGICAL KIT

271




Positioning Guide KIT Surgical Instruments

Guide Drill


- Initial drill : For initial drilling, assemble to the Single Guide to adjust the drilling depth
- Ø2.2 twist drill : Used with the bridge guide for initial drilling
- Ø3.0 twist drill : For subsequent drilling of Ø2.2 Twist Drill drilling path guide

D	Ø2.2	Ø3.0
Initial drill	GD2208NC	-
Twist drill	GD2213FNC	2D3011LC01



Single Guide

- Transparent material applied to facilitate the viewing of the position and direction for drilling
- 6 types with different mesiodistal crown diameters (Ø5.0~10.0)
- Packing unit : 2ea
- ※ Drilling depth adjusted to 6, 7 or 8mm using the marking line of the Initial Drill, based on the top line of the single guide
- ※ Disposable, Do not reuse




F5.0	F6.0	F7.0	F8.0	F9.0	F10.0
SGB050	SGB060	SGB070	SGB080	SGB090	SGB100

Guide Pin (Implant)

- Pin for checking the path and fixing the single guide in place after implant placement
- C = Connection

C	Mini	Regular
	PGSPM	PGSPR



Guide Pin

- A pin for checking the path around drilling and fixing the Single Guide in place



Bridge Guide

- Guide for adjusting the direction and distance for drilling
- Fan type : Selectable in 0.5mm increments (7~12.5mm)
- Compass type : Adjustable in 1 mm increments (5~24mm)
- Used after adjusting the distance in the distance setup of the mid panel of KIT



Type \ Distance	7~8.5	9~10.5	11~12.5	5~24
Fan	PGBRA070	PGBRA090	PGBRA110	-
Compass	-	-	-	PGBPA


Multi Joint Handle <sup>Option</sup>

- Instrument to place the guide from the outside of the mouth by connecting to the ball head of the Bridge Guide



Denture Guide <sup>Option</sup>

- Guide with adjustable angles for respective patients in edentulous cases
- Drilling is performed in the mouth with the angle fixed using an L-wrench after adjusting the angle according to the arch shape of the patient in a working model
- Marking line refers to the No. 2,3,4,5,6 positions of the teeth from the center



Positioning Guide KIT Surgical Instruments

L-wrench Option

- Instrument to adjust the size of the Denture Guide and fix it in place



LWC20H

Distance Setup Pin Option

- A pin for Bridge Guide compass type and denture guide fixation



PGDSP

SmartGuide KIT (OSGK)

Lower panel components

Guide Pin (4ea)  
SGP22



Round bur (2ea)  
RAHM1018




Twist Drill (2ea) (Ø2.2) For working models  
2D2208LC01



SGCB30S

X Sleeve  
(4ea)

Twist Drill  
(Ø2.2)  
SGTD2207S







Twist Drill  
(Ø3.0)  
QGTD3008



# SmartGuide KIT   Surgical Instruments

## SmartGuide

- Thermoplastic surgical guide
  - Freely deformable after immersion in about 70°C water for about 1 minute
  - Curing at room temperature after 1 minute from deformation
- ※ Disposable, Do not reuse; Use after low temperature disinfection (Do not autoclave or use hydrogen peroxide)

Type	Single	Free-end Bridge	2-Unit Br.: small	2-Unit Br.: large
				
	SGTSS	SGTFB90LS	SGTB63SS	SGTB85LS

## Twist Drill (Ø2.2) For working models

- Used for initial marking on the working model
- Number of use cycles: 10 times
- Additional drilling after using the round bur
- Recommended drilling speed: 1,200~1,500rpm

D	Ø2.2
	2D2208LC01



## Twist Drill

- A drill used through the guide in the mouth
- Enables stable drilling by connecting to the sleeve of SmartGuide
- After initial drilling with Ø2.2 drill, additional drilling with Ø3.0 drill is performed
- Recommended drilling speed : 1,200~1,500rpm

D	Ø2.2	Ø3.0
	SGTD2207S	QGTD3008



## Guide Pin

- Assembled to the working model for fixing the SmartGuide in place
- Connected to the SmartGuide sleeve

	SGP22
--	-------



## X Sleeve


- Instrument to check if the guide is produced as intended through CT scans or x-ray images by connecting to the SmartGuide sleeve
- After connecting to the SmartGuide outside the mouth, mount inside the oral cavity

	SGCB30S
--	---------



ESSET KIT (HESEK)

Applicable Products TSII / III SSII / III




• For ordering codes of single items of ESSET KIT, see pages 279-281

ESSET KIT Surgical Instruments

Crest Remover

- Marking the implant placement position after removing the narrow alveolar ridge horizontally
- Recommended drilling speed
  - Angled type : 1,200~1,500rpm
  - Straight type : 15,000~30,000rpm


L \ D	Ø5.0	Ø7.0
29	CERM50A	CERM70A
45	CERM50S	-



Twist Drill

- Marking the implant placement position
- Depth adjusted by assembling a stopper according to the implant length
- Recommended drilling speed : 1,200~1,500rpm


L \ TL \ D	Ø1.8
8.5 \ 33	2D1808LC01
10 \ 34.5	2D1810LC01
11 \ 36	2D1811LC01



Saw


- Incision of the narrowed ridge
- After vertical incision, perform the incision of the entire area in the mesial to distal directions
- Recommended drilling speed: 1,200~1,500rpm
- Recommended number of use cycles: 10 times
- Used after connecting to the saw protector at the position of the saw connecting groove

D	Ø7.0	Ø10.0	Ø13.0
	HSAW070	HSAW100	HSAW130



Expansion Drill

- Expansion of the ridge after incision
- Used in sequence according to the implant diameter  
F4.0 : I → II → III / F4.5 : I → II → III → IV
- Recommended drilling speed: 25~35rpm



Torque Wrench (Bar Type)

- Used for adjustment of placement position of implants and tightening of abutment and screws
- Torque is applied by pulling the bar and aligned to the line indicated with the torque value to be applied



TQWCB

Depth Gauge


- Instrument to release excessive torque by rotating the hex of the Expansion Drill with an open wrench when the hand piece does not move with the Expansion Drill stuck in alveolar bone in the process of removing the drill



ODG

Mount Extension

- Used for applying torque in manual mode in the process of placing or removing an Expansion Drill into the alveolar bone





ASMEL

EXP Mount Driver

- Used for applying torque for engine in the process of placing or removing an Expansion Drill into the alveolar bone

L	Type
Short (L)	AESMDS
Long (L)	AESMDL





• For ordering codes of single items of IM-Cure KIT, see pages 283-285

Metal Probe

- Instrument to measure the depth of periodontal disease
- Measuring periodontal pockets and identifying the shape of the periodontal pockets such as depth/size
- Marking line for probing in 1 mm increments



Plastic Probe

- Instrument to measure the depth of infection or periodontal disease around the implant
- Scratching of implant is prevented by using plastic material
- Flexible probe suitable for the curved form of alveolar bone
- Autoclave can be used
- Marking line for probing in 1 mm increments



Curette

- Instrument for removing subgingival sediments firmly attached to the granulation tissue of a specific area
- Gracey curette
- 01-02 : For removal of granulation tissue from anterior region
- 11-12 : For removal of granulation tissue from the mesial surface in anterior region
- 13-14 : For removal of granulation tissue from the distal surface in anterior region

Type	01-02	11-12	13-14
	HICC0102	HICC1112	HICC1314



IM-Cure KIT Surgical Instruments

Protect Screw


- Preventing infiltration of foreign substances into the internal connection of the implant using SmartBrush 2
- Tightened with 1.2 hex driver at force of 5Ncm



Type	Mini	Regular
	OPISB2PM	OPISB2PR

SmartBrush 1

- Used for peri-implantitis cleaning
- Used after connecting the Protect Screw to the implant after removing the patients prosthesis or abutment
- Recommended drilling speed: 1,200~1,500 rpm
- Recommended number of use cycles : About 1 minute per thread
  - ※ Do not use for longer than 4 minutes
- Be sure to polish with saline irrigation and suction
  - ※ Disposable and do not reuse (Must be discarded after use)



L	
Short	OPISB1S
Long	OPISB1L

SmartBrush 2

- Used for peri-implantitis cleaning
- Used after connecting the Protect Screw to the implant after removing the patients prosthesis or abutment
- Be sure to polish with saline irrigation
- Recommended drilling speed: 1,200~1,500rpm
- Recommended number of use cycles : 1~2 minutes
  - ※ Excessive use for longer than 3 minutes may result in fracture or bending of the product
  - ※ Disposable, Do not reuse (Must be discarded after use)



L \ D	F3.0 / F3.5	F4.0 / F4.5	F5.0 / F5.5	F6.0	F7.0
Short	OPISB23S	OPISB24S	OPISB25S	OPISB26S	OPISB27S
Long	OPISB23L	OPISB24L	OPISB25L	OPISB26L	OPISB27L

Metal Scaler

- Used for removing plaque or foreign substance from the surface of the implant by connecting to an ultrasonic scaler
- Used as a secondary instrument after using SmartBrush 1 or SmartBrush 2
- Bendable tip of the product for easy access
- EMS, KaVo and SATELEC types available

Type	EMS	KaVo	SATELEC
	HICSME	HICSMK	HICSMS



Plastic Scaler Connector

- Used by assembling to a plastic scaler tip
- Do not use for removing foreign substances from the implant surface
- EMS, KaVo and SATELEC types available
- A = Angle


A \ Type	EMS	KaVo	SATELEC
125°	HICSCE	HICSCK	HICSCS
100°	HICSCE2	HICSCK2	HICSCS2



Plastic Scaler Tip

- Used for removing foreign substances from the abutment or crown by connecting to a SmartScaler
  - ※ Do not use on the implant surface
- Packing unit : 30ea/1set

PSTP
------



ESR KIT Easy Screw Removal KIT (OESRK)


ESR Full KIT Easy Screw Removal Full KIT (OESRFK)

• Including the same components as ESR KIT and allows holding of the components provided by other companies

Applicable Products Nobel Biocare Active/Replace / Straumann Bone Level / Astra Osseo Speed TX  
3i Full OSSEOTITE Tapered Certain / Zimmer Tapered / Biohorizons Internal


SURGICAL KIT

286



SURGICAL KIT

287





Items not included in the KIT

Guide								
Nobel	Active	Replace	3i	Tapered Certain		Straumann	Bone Level	Roxolid SLActive
	OGNA01L	OGNR02L		OGIF01L	OGS01L			
	OGNA02L	OGNR03L		OGIF02L	OGS02L			
		OGNR04L						OGSTRS
							OGS02L	OGSTRL
Astra	Osseo Speed TX		Biohorizons	Internal	External	Zimmer	Tapered	
	OGAO01L			OGZB01L	OGBES		OGZB01L	
	OGAO02L			OGZB02L	OGBEL		OGZB02L	
SR Drill			SR Tip		Screw Holder		Re-Tap	
OSRD09L			ORT18L		OSHR18L		ORTR18L	



Guide

- Used for centering and prevention of shaking of SR Drill, SR Tip, etc. by connecting and fixing to the implant
- Use according to implant type and diameter  
(Internal/submerged type products of 6 overseas manufacturers)
- Short or Long types selected according to the intermaxillary distance
- ■ Used in common
- C = Connection / D = Diameter


Osstem

C \ Type	TS		SS	
	Short	Long	Short	Long
				
Mini	OGTMS	OGTML	OGUMS	OGUML
Regular	OGTRS	OGTRL	OGSRS	OGSRL
Wide	-	-	OGSRS	OGSRL


Nobel Biocare


D \ Type	Active		Replace	
	Short	Long	Short	Long
				
ø3.5	OGNA01S	OGNA01L	-	-
ø4.3	OGNA02S	OGNA02L	OGNR02S	OGNR02L
ø5.0	OGNA02S	OGNA02L	OGNR03S	OGNR03L
ø6.0	-	-	OGNR04S	OGNR04L

Nobel Biocare


D \ Type	MkIII	
	Short	Long
		
ø3.3	OGUMS	OGUML
ø3.75	OGURS	OGURL
ø4.0	OGURS	OGURL
ø5.0	OGUWS	OGUWL

Straumann



D \ Type	Bone Level	
	Short	Long
		
NC (3.3)	OGSB01S	OGSB01L
RC (4.1)	OGSB02S	OGSB02L
RC (4.8)	OGSB02S	OGSB02L

D \ Type	Roxolid SLActive	
	Short	Long
		
RN (3.3)	OGSTRS	OGSTRL
RN (4.1)	OGSTRS	OGSTRL
RN (4.8)	OGSTRS	OGSTRL
WN (4.8)	OGSTRS	OGSTRL

Astra

D \ Type	Osseo Speed TX	
	Short	Long
		
Small (3.5 s)	OGAO01S	OGAO01L
Small (4.0 s)	OGAO01S	OGAO01L
Large (4.5)	OGAO02S	OGAO02L
Large (5.0)	OGAO02S	OGAO02L
Large (5.0 s)	OGAO02S	OGAO02L


3i

D \ Type	Full Osseotite Tapered Certain		D \ Type	Full Osseotite Tapered	
	Short	Long		Short	Long
					
3.25	OGIF01S	OGIF01L	Ø4.0	OGURS	OGURL
4.0	OGIF02S	OGIF02L	Ø5.0	OGURS	OGURL
5.0	OGIF02S	OGIF02L	Ø6.0	OGURS	OGURL
6.0	OGIF02S	OGIF02L			


Reverse Drill

- Instrument used for removing fractured screws
  - Be sure to use with a suitable guide for the implant
  - When the red marking of the reverse driver is shown above the guide assembled to the implant, use a screw holder to remove the fractured screw
  - For hand mode / Rotating direction : Reverse rotation / Number of use cycles: 10 times
- ※ Do not use more than 10 times; Do not reuse

L \ Type	M1.6	M1.8	M2.0
Short	-	ARVDRSC	ARVDRSC
Long	ARVDMLC	ARVDRLC	ARVDRLC




Zimmer

D \ Type	Tapered	
	Short	Long
		
Green (3.7)	OGZB01S	OGZB01L
Green (4.1)	OGZB01S	OGZB01L
Green (4.7)	OGZB01S	OGZB01L
Green (6.0)	OGZB02S	OGZB02L



Screw Removal Drill (SR Drill)

- Used for removal to form a hole in fractured screws
  - Be sure to assemble to the guide and remove the cut chips by suction with irrigation to the window
  - Select Short or Long types according to the intermaxillary distance
  - Drilling until the red line around the handle is not visible
  - Recommended drilling speed: 1,200~1,500rpm in reverse rotation / Number of use cycles: 5 times
- ※ Be sure to use with a guide assembled / Do not exert excessive vertical force / Do not soak in hydrogen peroxide
- ※ Disposable, Do not reuse
- Short : Sold as a single item

L \ Type	M1.6	M1.8	M2.0
Short	OSRD08S	OSRD09S	OSRD10S
Long	OSRD08L	OSRD09L	OSRD10L



Biohorizons

D \ Type	Internal (Tapered Bone Level)		D \ Type	External	
	Short	Long		Short	Long
					
Yellow	OGZB01S	OGZB01L	Ø3.5	OGUMS	OGUML
Green	OGZB01S	OGZB01L	Ø4.0	OGURS	OGURL
Blue	OGZB02S	OGZB02L	Ø5.0	OGBES	OGBEL
			Ø6.0	OGBES	OGBEL

Torque Driver Handle


- Used by rotating by hand after assembling with products such as SR tip, AR tip, and screw holder



Reverse Driver

- Instrument used for removing fractured screws
  - Be sure to use with a suitable guide for the implant
  - When the red marking of the reverse driver is shown above the guide assembled to the implant, use a screw holder to remove the fractured screw
  - For hand mode / Rotating direction : Reverse rotation / Number of use cycles : 10 times
- ※ Do not use more than 10 times
- C = Connection


L \ C	Mini	Regular/Wide
Short	-	ORVDRS
Long	ORVDML	ORVDRL



Re-tap

- Instrument to restore the thread to the initial state when the screws cannot be fastened due to damage to the internal thread of the implant
- Thread formed in hand mode with a Torque Wrench or Ratchet Wrench

Type	M1.6	M1.8	M2.0
	ORTM	ORTR18	ORTR



Screw Removal Tip (SR Tip)

- Used for removing fractured screws by rotating the screw removal tip in the hole on the fractured surface of the screws formed by using the screw removal drill (SR Drill)
  - Rotating direction : Reverse rotation
- ※ Disposable, Do not reuse

L \ Type	M1.6	M1.8	M2.0
Short	ORT16S	ORT18S	ORTS
Long	ORT16L	ORT18L	ORTL



ESR Handle

- Instrument to fix the guide to the implant


Type	OARH
------	------



Screw Holder

- Removing partially protruding fractured screws by assembling with a screw holder
- Color coded for easily visible indication of types
- Rotating direction : Reverse rotation


Type	M1.6	M1.8	M2.0
	OSHM	OSHR18	OSHR



Abutment Removal Tip (AR Tip)

- Used for fractured abutment, mount partially remaining and stuck in the implant
- The AR Tip is assembled with the fractured abutment hole, fixed in place tightly and abutment is removed using forceps, etc.
- Mini : removing screws with a slipped hex
- The Mini AR Tip is assembled to the slipped hex, and rotated in the reverse direction to connect to the screw for removal

L \ Type	Mini	Regular
Short	OARTMS	OARTRS
Long	OARTML	OARTRL
Ex.Long	OARTMEL	OARTREL



Slot Driver

- Instrument to be used by forming a slot with Ø0.8 bur, when force cannot be exerted using a driver due to the damaged hex of Healing Abutment, Cover Screw, or Abutment Screw.



OTSD07

Transfer Abutment Separate Tool

- Used for releasing the jamming caused by Non-hex Transfer Abutment stuck from the contact of the implant and the morse taper
- Commonly used for both Mini and Regular: the body end is used for Mini and for Regular, it is placed into the 2-stage groove
- For ease of separation, the separate tool body is placed into the inner hole of the abutment after removing the abutment screw, and the driver is rotated forward to integrate the body and the abutment.
- If separation is difficult, use the tool after connecting a Ratchet Wrench to the driver



Driver



Body

Driver

TASD

Body

TASB

Set

TAST

OSSTEM<sup>®</sup>  
IMPLANT

EFR KIT Easy Fixture Removal KIT (OSFRK)




Top panel components

Implant Wrench  
FRDFE



Applicable Products **TSII / III** **SSII / III** **Ultra-wide**



EFR Full KIT Easy Fixture Removal Full KIT (OSFRFK)

• Including the same components as EFR KIT, which can hold the components provided by other companies

Top panel components

Implant Wrench  
FRDFE



Torque Wrench  
TW400B



Applicable Products

**Nobel Biocare** Active/Replace / **Straumann** Bone Level / **Astra** Osseo Speed TX  
**3i** Full OSSEOTITE Tapered Certain / **Zimmer** Tapered / **Biohorizons** Internal



EFR Full KIT    Surgical Instruments

Items not included in the KIT

Remover Screw								
Nobel	Active		Replace					
	Normal	Fracture	Normal	Fracture				
	FRSMNA35	FRSMNA35F	FRSMNR35	FRSMNR35F				
	FRSR40	FRSR40F	FRSR40	FRSR40F				
	FRSW50	FRSW50F	FRSW50	FRSW50F				
Straumann	Bone Level		3i	Full Osseotite Tapered Certain		Biohorizons	Internal	
	Normal	Fracture		Normal	Fracture		Normal	Fracture
	FRSM33	FRSM33F		FRSMI325	FRSMI325F		FRSRZ41	FRSRZ41F
	FRSRS41	FRSRS41F		FRSRI40	FRSRI40F		FRSWZ47	FRSWB46F
	FRSWS48	FRSWS48F		FRSWI50	FRSWI50F		FRSWZ60	FRSWB46F
Zimmer	Tapered		Astra	Osseo Speed TX		Remover Body		
	Normal	Fracture		Normal	Fracture	FRBW57S		
	FRSMZ37	FRSMZ37F		FRSMNA35	FRSMNA35F	FRBW57L		
	FRSRZ41	FRSRZ41F		FRSRA40	FRSRA40F	FRBUW60S		
	FRSWZ47	FRSWZ47F		FRSR40	FRSR40F	FRBUW60L		
	FRSWZ60	FRSWZ47F		FRSW50	FRSW50F			

Remover Screw

- Acting as a support structure for reverse rotation of the remover body after connected and fixed to the implant
- Used according to the type and diameter of the implant to be removed (Internal/submerged type products of 6 overseas companies, Normal/Fracture)
- Fracture mode is used for removing implants with the hex completely fractured
- Compatible with products of 6 overseas companies
- Recommended tightening torque: Regular/Wide 80Ncm, Mini 60Ncm
- T = Type    ※ Disposable, Do not reuse



Osstem

T \	Mode	Mini Ø3.5 /-	Regular Ø4.0~4.5 / P4.8	Wide Ø5.0 / P6.0
TS/SS	Normal	FRSM35	FRSR40	FRSW50
	Fracture	FRSM35F	FRSR40F	FRSW50F

Nobel Biocare

T \	Mode	Mini Ø3.5	Regular Ø4.3	Wide Ø5.0/6.0
Active	Normal	FRSMNA35	FRSR40	FRSW50
	Fracture	FRSMNA35F	FRSR40F	FRSW50F
Replace	Normal	FRSMNR35	FRSR40	FRSW50
	Fracture	FRSMNR35F	FRSR40F	FRSW50F

Straumann

T \	Mode	Mini Ø3.3	Regular Ø4.1	Wide Ø4.8
Bone Level	Normal	FRSMS33	FRSRS41	FRSWS48
	Fracture	FRSMS33F	FRSRS41F	FRSWS48F

Astra

T \	Mode	Mini Ø3.5	Regular Ø4.0	Regular Ø4.5	Wide Ø5.0
Osseo Speed TX	Normal	FRSMNA35	FRSRA40	FRSR40	FRSW50
	Fracture	FRSMNA35F	FRSRA40F	FRSR40F	FRSW50F

3i

T \	Mode	Mini Ø3.25	Regular Ø4.0	Wide Ø5.0/6.0
Full Osseotite Tapered Certain	Normal	FRSMI325	FRSRI40	FRSWI50
	Fracture	FRSMI325F	FRSRI40F	FRSWI50F

Zimmer

T \	Mode	Mini Ø3.7	Regular Ø4.1	Wide Ø4.7	Ultra-wide Ø6.0
Tapered	Normal	FRSMZ37	FRSRZ41	FRSWZ47	FRSWZ60
	Fracture	FRSMZ37F	FRSRZ41F	FRSWZ47F	FRSWZ47F

Biohorizons

T \	Mode	Mini Ø3.8	Regular Ø4.6	Wide Ø5.8
Internal	Normal	FRSRZ41	FRSWZ47	FRSWZ60
	Fracture	FRSRZ41F	FRSWB46F	FRSWB46F

EFR Full KIT Surgical Instruments

Screw Driver

- Driver to connect and fix the remover screw to the implant
- Remover screw Recommended tightening torque: Regular/Wide 80Ncm, Mini 60Ncm
- The same type as the remover screw is selected
- T=Type

T	Mini	Regular	Wide
	FRSDM23	FRSDR25	FRSDW30



Torque Wrench

- Used to remove the implant with the remover body after tightening with a screw driver
- Torque applied up to 400Ncm (60/80/200/300/400Ncm scale display)
- Torque applied by aligning the center of the bar with the torque value to be applied by pulling the bar
- Washed and sterilized after use for storing

	TW400B
--	--------



Remover Body

- Instrument to exert torque in the implant loosening direction by connecting to a remover screw
- Used according to the diameter of the implant to be removed
- ※ Disposable, Do not reuse
- The same type as the remover screw is selected
- T=Type

	T	Mini	Regular	Osstem only Wide	For overseas other companies Wide	Ultra-wide
Short		FRBM35S	FRBR40S	FRBW50S	FRBW57S	FRBUW60S
Long		FRBM35L	FRBR40L	FRBW50L	FRBW57L	FRBUW60L



Implant Wrench

- Wrench to remove implant from the remover body

	FRDFE
--	-------



Torque Extension


- Extension of the length of screw driver and remover body (up to 10mm)

	OTE
--	-----



# Dr. Cho's Instrument KIT (DCHOKIT)

- Optimal implant surgery KIT based on years of clinical know-hows
- Composed of 10 types of instruments (1ea each)



# Osstem Basic Instrument KIT (OBKIT)

- Commonly used implant surgery KIT
- Composed of 25 types of instruments (1ea each)



Osstem Basic Instrument KIT (OBKIT)

SURGICAL KIT

304

**Periosteal Elevator (Selden)**

- Retracted and fixed for gingival tissue flap
- W : 10/13mm

EP23-W

**Caliper**

- Castroviejo

LPC90-W

**Bone Well**

BWSUS1-W

**Scissor (LaGrange)**

- Compound (curved)
- L : 115mm (±5)

SCLC115-W

**Surgical Curettes (Gracey)**

CGR11-12-W

**Surgical Curettes (Surgical Curettes, CM11)**

URCM11-W

**Minesota Retractor**

RTCRM-W

**Towel Clamp**

- Towel Clamp, Backhaus
- L : 135mm (±5)

CPTC135-W

SURGICAL KIT

305

**Periosteal Elevator (MOLT9)**

- Lifting mucosal periosteum after gingival tissue incision
- W : 8.2/4.2mm

EP9-W

**Tweezer (Wide)**

- L : 155mm

PCW150-W

**Bone Rongeurs**

- Friedman
- L : 140mm (±5)

RNGF140-W

**Surgical Curettes (Surgical Curettes, CM10)**

URCM10-W

**Periosteal elevator (Prichard)**

- W : 11/4.9mm

EP9R3-W

**Scalpel Handie (Straight Type)**

SHS-W

**Scissor (Tissue Scissor)**

- Straight
- L : 150mm (±5)

SCTC115-W


**Mallet**

- Autoclave can be used

ML25-W


# Custom KIT (OCTK)

- KIT used to disinfect some of the surgical instruments or to store new spare tools
- Additional 3 types of rubber (large, medium, small), which can be used according to user preference
- Sterilizable material (132℃, 15 minutes)




# Healing Case (OHAC)

- A case for temporary storage and cleaning of Healing Abutment during the prosthesis procedure
    - Upper prosthesis for additional mounting: transfer / temporary / angled / cover screw / pick-up & transfer impression coping/ OB anchor/ temporary crown (Only the Healing Abutment can be assembled with the top plate)
  - Similar to the tooth arrangement, a total of 28 cells are arranged with 7 cells each in the upper / lower and left / right sections
  - Sterilizable material (132℃, 15 minutes); sterilization required for reusing the case
- ※ This product is not a case for reuse of Healing Abutment




# Osteo KIT (OSTK)

- Concave osteotome : KIT used for sinus lift procedure(maxillary sinus floor elevation) to vertically increase the amount of alveolar bone available in the maxillary anterior region
- Expanding osteotome : KIT used to increase the primary stability of the implant in low quality bones by densifying the trabeculae of bone while preserving the bone instead of removing it
- Stopper for adjusting the depth of procedure




# Osteotome KIT (AOST)

- KIT used for sinus lift procedure(maxillary sinus floor elevation) to vertically increase the amount of alveolar bone available in the maxillary anterior region
- Included in concave type only
- Stopper for adjusting the depth of procedure





Sinus KIT (ASLK)

- KIT containing various tools for maxillary sinus floor elevation (sinus lift procedure)
- Lateral approach instrument for sinus procedure
- Components (5 types)
  - Freer elevator : OFE
  - Bone graft carrier : OBGC
  - Membrane separator (circle type) : OMSC
  - Sinus curette-short : OSCS
  - Sinus curette-long : OSCL




Bone Spreader KIT (OBSOK)

- KIT used for expanding narrowed alveolar ridge
- Offset type for easy operation
- Components (4 types)
  - OBSO22F, OBSO28F, OBSO35F, OBSO35R



- Use for alveolar bone expansion
- Offset type for easy operation
- Depth marking corresponding to the implant length






		(Unit : mm)			
Code	Spec.	Tip length			
		7	8.5	10	11.5
OBSO22F	Thickness	1.15	1.3	1.45	1.6
	Width	2.1	2.2	2.2	2.2
OBSO28F	Thickness	1.15	1.3	1.45	1.6
	Width	2.65	2.8	2.8	2.8
OBSO35F	Thickness	1.3	1.45	1.6	1.8
	Width	3.3	3.5	3.5	3.5
OBSO35R (round type)	Thickness	1.85	2.1	2.3	2.55
	Width	3.3	3.5	3.5	3.5

Ridge Split KIT Straight (ORSSK)

Straight

- Chisel : Used for expanding narrowed alveolar ridge
- Blade holder : Malletting (as seen below) enabled by tightening a #15 blade when it is difficult to make a bone incision using bur due to low bone quality
- Components
  - Ridge split chisel : ORSS15, ORSS20, ORSS25, ORSS30
  - Blade holder : ORSBH






		Tip length			
Code	Spec.	7	8.5	10	11.5
ORSS15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSS20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSS25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSS30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4

Ridge Split KIT Offset (ORSOK)

Offset

- Chisel : Used for expanding narrowed alveolar ridge
- Blade holder : Malletting enabled by tightening a #15 blade when it is difficult to make a bone incision using bur due to low bone quality
- Components
  - Ridge split chisel : ORSO15, ORSO20, ORSO25, ORSO30
  - Blade holder : ORSBH



		Tip length			
Code	Spec.	7	8.5	10	11.5
ORSO15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSO20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSO25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSO30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4

**OSSTEM<sup>®</sup>**  
IMPLANT