

ATLANTIS ISUS™ compatibility chart

– Implant systems

Extended Assortment



Manufacturer/ Company	Implant system	Implant level	Abutment level	Screw Order REF number from Manufacturer/Company*
Alphatech	INTERNAL			
	D3.4, D3.8, D4.3 Implant Level	x		188705 (5 screws)
	D5, D6, D7 Implant Level	x		188705 (5 screws)
	EXTERNAL			
Alpha Bio	Abutment level		x	188704 (5 screws)
	INTERNAL			
	Internal Hex	x		STLAS
	EXTERNAL			
BEGO	TSA Abutment		x	SF
	INTERNAL			
	Semados S/RI D3.25/D3.75	x		53999
	Semados S/RI D4.1, D4.5, D5.5	x		53999
Biotech	Semados Mini	x		56983
	INTERNAL			
	Bis D5/D5.4 Abutment Level		x	E007
	Bis D3/D3.4/D3.6/D3.9/D4.4 Abutment Level		x	E007
BT	Conical Abutment		x	N/A
	INTERNAL			
BTI	Lock D4.5	x		370 BTIVPT
	INTERNAL			
	Tiny D3.5	x		TTTH
	Interna Universal D4.1	x		INTTUH
Bredent	Interna Wide D5.5	x		INTTUH
	Multi-IM Abutment		x	TTMI
	INTERNAL			
	SKY D4.1 Implant Level	x		SKY PS22 (6X)
Dentaurum	SKY Fast/Fixed Abutment Level		x	SKYFFSPK (6X)
	INTERNAL			
	TioLogic S Implant Level	x		386-910-00
	TioLogic M Implant Level	x		386-910-00
Dentegris	TioLogic L Implant Level	x		386-910-00
	TioLogic Abutment Level		x	386-912-00
	Tiox D3.5, D4.5 Implant Level	x		330-790-00 (5mm)
	INTERNAL			
DIO	D3.3 Implant Level	x		ZS 01
	D3.75 Implant Level	x		ZS 01
	D4.5/D4.75 Implant Level	x		ZS 01
	D5.5 Implant Level	x		ZS 03
Dental Ratio	INTERNAL			
	OKTAGON Implant Level	x		84960
	INTERNAL			
	Legacy D3	x		1000-04
Implant Direct	Legacy D3.5 / 4.5 / 5.7	x		1000-01
	S-plant D3.7, D4.7, D5.7	x		1000-01
	S-InDirect D5.0	x		1000-01
	Overdenture Screw-Retained Abutment		x	1000-05
Imtec	INTERNAL			
	Implant-Level D3.5/D4.3	x		5052HX
	Implant-Level D5.1	x		5052HX
Medical Instinct	INTERNAL			
	BoneTrust D3.4 Implant Level	x		160-100004
	BoneTrust D4, D5 Implant Level	x		160-100004
	D3.4 /D4 Direct Abutment		x	165-344055

* **Note:** We recommend to check Screw REF number with the manufacturer.
All trademarks and company names are the property of their respective owners.

Extended Assortment

Manufacturer/ Company	Implant system	Implant level	Abutment level	Screw Order REF number from Manufacturer/Company*
Megagen	INTERNAL			
	Exfeel Internal Implant Level	x		EXIMS100
	Exfeel Internal Abutment Level		x	IRCS200
	Rescue Internal Implant Level	x		ARAS250L
	EZ Plus Internal Implant Level Small	x		EPAS100
	EZ Plus Internal Implant Level Regular/Wide	x		AMAS100
Medentis	INTERNAL			
	ICX Templant	x		C-007-000001
MIS	INTERNAL			
	Direct Standard	x		MD-S0220
	Direct Wide	x		MD-S0220
	Multiunit		x	MM-S0330
Neoss	INTERNAL			
	Narrow D3.25 Implant Level	x		31240
	D4.1 Implant Level	x		31150
	D4.1 Access Abutment		x	31229
Osseolink	INTERNAL			
	D3.5/D4 /D5 Implant Level	x		UNIPS
Osstem	INTERNAL			
	GS Mini	x		GSASM
	GS Standard	x		GSASR
	GS Convertible Mini Abutment D4.0		x	GSFSM
	GS Convertible Standard Abutment D4.8		x	GSFSR
	GS Convertible Wide Abutment D6.0		x	GSFSR
	US Mini	x		ASM200
Pitt Easy	INTERNAL			
	D3.25, D3.75, D4, D4.9, D6.5	x		90276
Platon	INTERNAL			
	PF D4.4 Implant Level	x		QPB-0001
	PF D4.8 Implant Level	x		QPB-W001
	Pro PF D4.1, D5.0 Implant Level	x		QPB-0001
	PF D4.1/D4.4 Abutment Level		x	QPA-0005/QPB-0001
	PF D4.8 Abutment Level		x	QPA-W005/QPB-W001
POI	INTERNAL			
	EX D3.7 Implant Level	x		CDO12600
	EX D4.2 Implant Level	x		CDO12650
	EX D.7 Implant Level	x		CDO12700
	EX D5.2 Implant Level	x		CDO12750
Prowital	INTERNAL			
	Wieland D3.5, D4.3, D5 Implant Level	x		1.213.1600
Renew Biocare	INTERNAL			
	UMA Abutment		x	11-000052
SIC	INTERNAL			
	D3.3, D4.2, D5.3	x		936658
Simpler	INTERNAL			
	Bar Abutment		x	AS 21001
Serf	INTERNAL			
	EVL Plus D3.3	x		VPEG 3,3-DP30700050
	EVL Plus D4, D5	x		VPEG 4 ET 5-DP30700070
Schuetz	INTERNAL			
	Dental Impla D3.3, D4.2, D5.3	x		636649
Sweden Martina	INTERNAL			
	Global D3.8, D4.3, D4.8, D5.5	x		U-VMOU1.8
	PAD Abutment Level		x	PAD-VP-140
	Premium/Kohno D3.3, D3.8	x		VM2-180
	Premium/Kohno D4.25, D5	x		VM2-200
Tekka	INTERNAL			
	Progress D4.5 Implant Level	x		DVPQP
	Progress D5, D6 Implant Level	x		DVPQP
	Krestal D3.6 Implant Level	x		DVPQP
	Monobloc D3.4	x		DVPIP
	Monobloc D4, D4.5, D5	x		DVPIP
	Quattro	x		DVPAOFF
Thommen Medical	INTERNAL			
	SPI D3.5 Implant Level	x		4.03.500
	SPI D4, D4.5, D5, D6 Implant Level	x		4.03.501
	SPI Abutment Level		x	4.03.506
	SPI VARIO multi Abutment		x	4.03.507

* **Note:** We recommend to check Screw REF number with the manufacturer.
All trademarks and company names are the property of their respective owners.

ATLANTIS ISUS™ Turnaround times

ATLANTIS ISUS™ implant suprastructures are available for major implant systems, in titanium and cobalt-chrome. Turnaround times vary depending on material.

When the case is received at the ATLANTIS ISUS™ production facility and the manufacturing process has started, your ATLANTIS™ WebOrder account will be updated with the case estimated shipment date. Standard turnaround times are as follows:

ATLANTIS ISUS™ in-house production time (days)		
	Titanium	Cobalt-Chrome
ATLANTIS ISUS™ bridge	5	7
ATLANTIS ISUS™ hybrid	5	7
ATLANTIS ISUS™ bar	5	7
ATLANTIS ISUS™ 2 in1	12	N/A

Determine return date

1. Select your choice of delivery method and calculate the number of days it will take to ship the case from the dental laboratory to the ATLANTIS ISUS™ production facility.
2. Add the days the case will be in-house at the ATLANTIS ISUS™ production facility. Calculation of in-house production time starts when the case is received at the ATLANTIS ISUS™ production facility before noon.
3. DENTSPLY Implants will ship the case free of charge via next day delivery.

Note:

- Design approvals to process case must be granted by 6 pm CET.
- Material, on-hold and case information issues must be resolved by midnight CET.
- Weekends and holidays are excluded.
- Time needed for **Design approval** or resolving **On hold** information may delay the case.

Example

Order and all case materials for an ATLANTIS ISUS™ bridge in titanium is received on Monday before noon. Design for approval will be sent out.

If approval takes place the same day by 6 pm CET, the case will be shipped to the dental laboratory on Friday.

Additional information

The **ATLANTIS ISUS™ prosthetic screw** is offered for selected major implant systems. The screws can be ordered together with the ATLANTIS ISUS™ implant suprastructure in the ATLANTIS™ WebOrder and with the ATLANTIS ISUS™ prescription form.

The ATLANTIS ISUS™ implant suprastructures are always compatible with the implant/abutment screw from respective manufacturers system. If the ATLANTIS ISUS™ prosthetic screw is not available, the screws can be ordered from the implant system manufacturer. As a service all screws are listed with the manufacturers REF number in the ATLANTIS ISUS™ compatibility chart. Note that we recommend to check the REF number with the implant system manufacturer.

For the most current version of the ATLANTIS ISUS™ implant compatibility chart and turnaround times, visit www.dentsplyimplants.com or your ATLANTIS™ WebOrder account.

